



airwalk se

original directions for use

Operating instruction and service manual
Body weight support / unweighting device
h/p/cosmos airwalk se
manual version 1.3

Development, production, sales & service

h/p/cosmos sports & medical gmbh
Am Sportplatz 8
DE 83365 Nussdorf-Traunstein
Germany
phone +49 86 69 86 42 0
fax +49 86 69 86 42 49
email@h-p-cosmos.com
www.h-p-cosmos.com

© 2014 h/p/cosmos sports & medical gmbh
All rights reserved

This manual is only valid for the firmware-/software version noted on the first page of this manual and only for the original configuration of the first delivery of the machine.

Firmwareupdates, softwareupdates, changes of the system configuration or retrofittings of additional equipment or accessories can result in invalidity of this manual. In case of alterations of the device or the additional equipment always the actual version of the manual or the according additional information has to be considered.

[1.] Introduction

Dear customer,

We would like to express our gratitude for putting your trust in us, in deciding for this high standard the unweighting system. Since 1988 we have been developing and manufacturing running machines for sports and medical application. When it comes down to technology, design and safety, we have set extremely high standards for ourselves.

Because the unweighting device is a motor-driven and powerful device, you must pay special attention to the mentioned safety regulations. If proper notice is taken of the safety regulations, the operation of our unweighting device is almost without any risk. The neglect of the safety regulations could result in dangerous situations and serious accidents. Therefore please carefully read the installation and operation manual and the danger precautions before taking the device into operation.

Simple maintenance, if described, can easily be done by yourself. We recommend to call our competent service team or entering into a maintenance contract for a routine service in an interval of 6 or 12 months. A form for registration of your institution and device is included in the delivery folder. In order to be able to supply you with the latest technical information and service, it is important for you to fill out the form. Please fill out the form for registration immediately and send it back via fax.

The instruction manual as a firm part of the delivery has to be accessible for all users at any time. This instruction manual has been written with great care. Should you, however, still find any details, which do not correspond with your device, please give us notice, so that we can correct any mistakes as soon as possible. Subject to alterations without prior notice,

We wish you a lot of fun and success while exercising and working with your h/p/cosmos unweighting device.




Franz Harrer

President

h/p/cosmos sports & medical gmbh

Comment:

The unweighting device h/p/cosmos airwalk se is compatible with below listed running-machine series.



h/p/cosmos 150/50 locomotion DE med



h/p/cosmos mercury



h/p/cosmos quasar



h/p/cosmos pulsar



h/p/cosmos venus 200/75if mounted in a pit

[2.] Content

[1.] Introduction	3
[2.] Content	5
[3rd] Device description	9
[4.] Safety, warnings, precautions, prohibitions	10
[3.A] Safety precautions, safety regulations, prohibitions and warnings	10
[3.B] Danger zones	10
[3.C] Emergency release	11
[3.D] General Instructions / Safety	11
[3.E] Instructions for safety and operation	13
[3.E1] Electric safety / Safety classes	13
[3.E2] Patient security according to VDE 0750 (IEC 601-1) type B	13
[3.E3] Potential equalization	13
[3.E4] Connection of units installed outside medically utilized locations	14
[3.E5] Medical electrical systems	14
[3.E6] Subject's surroundings	14
[3.E7] Protection against the danger of electrical shocks	15
[3.E8] Environmental requirements	15
[3.F] Economic life-time	16
[3.G] Symbols used	16
[3.H] Safety standards, norms	16
[3.H1] The CE mark	17
[3.H2] EMC Electromagnetic Compatibility	17
[3.H3] General advice	17
[3.I] Fields of application / intended use	18
[3.I1] Professional application in sports and fitness	18
[3.I2] Professional application in the medical field	18
[3.I3] h/p/cosmos devices - MDD medical device directive. Safety standards	18
[3.I4] Forbidden use - contraindications	21
[5.] Installation & commissioning	23
[5.A] Transport, unpacking and packaging	23
[5thB] Transport to upper floors and through narrow doors	23
[5.C] Mechanical installation	24
[5.D] Electrical installation	25
[5.D1] Name plate	26
[5thD2] Electric safety measurements and „First Measured Values“	26
[5.D3] Potential equalization	26
[5.E] Starting the unweighting device	28
[5.E1] Switching the device on	28
[5.E2] Switching the device off	28
[5.F] Installation, Overview, Checklist	28

[6.] Operation	37
[6thA] Controlements and displays: Manual control unit and housing unweighting unit	37
[6thA1] Manual control unit of Levi unweighting unit	37
[6thA2] Displays of the Levi unweighting unit and the connection for manual control	38
[7.] Training	39
[7.A] General notes	39
[7.A1] Suitable clothing	39
[7.A2] Verwendete Laufbänder und Rollstuhlrampen:	39
[7.A3] Before starting training	39
[7.B] h/p/cosmos airwalk Vests	40
[7.B1] Putting on the Vest	40
[7.B2] Sizes	40
[7.B3] Adjust Chest Straps	41
[7.B4] Buckle Waist	42
[7.B5] Cuff Leg	43
[7.B6] Adjust Leg Straps	43
[7thC] Attaching the Vest to the Unweighting Bar	44
[7.C1] Correcting Vest Problems	44
[7.D] Velcro-Gait Control Straps	44
[7thD1] Internal Foot Rotation, Excessive Stride Length	45
[7.D2] External Foot rotation	45
[7thE] Putting on the h/p/cosmos airwalk Vest for Wheelchair Patients	46
[7thF] Entering the unweighting device – supportign weight / standing up	47
[7thG] Setting the support weight and range of motion	47
[7.H] General advice for training	49
[7.I] Unweighting / weight support	49
[7.J] speed	49
[7.K] Elevation	50
[7.L] Trainingduration/ distance und training frequency	50
[7.M] Training troubleshooting	50
[7.N] After training	51
[7thO] At the end of a training day	51
[7.P] Patient documentation	51
[7.Q] Patient becomes unconscious	51
[8.] Maintenance and safety inspections	52
[8.A] Preventiv maintenance	53
[8.B] Immediate maintenance	53
[8.C] Regular inspections/examinations	53
[8thC1] Visual inspection dirt / damage – daily before training	54
[8thC2] Emergency release unweighting device – every 3 months	54
[8thC3] Running machine – follow manufacturers guidelines	54
[8thC4] Changing the rope – min. every year	55
[8thC5] Spare parts	56
[8thC6] Protective Earth Resistance RPE measurement	56

[8.C7] Isolation Resistance R_{ISO} measurement	57
[8thC8] Equivalent (alternative) leakage current I_{EDL} measurement	57
[8thC9] Electric Safety Tester an measurements	57
[8thC10] Building installation: Electric Checks, Protective Earth Function, RCD, Leakage Current	58
[9.] Trouble shooting.....	59
[9.A] Mechanical / Noise problems	59
[9.B] Interference factor	59
[9.B1] Electrostatic Discharge	59
[9.B2] Source of Interference	60
[9thC] Voltage at the device housing / electric shock	60
[9thC1] Open (interrupted) earth wire / ground – shown here on the esample of a running machine	60
[9.D] Electrical interference	61
[10.] Technical specifications	62
[10thA] h/p/cosmos airwalk se	62
[10thB] EMC electromagnetic compatibility. Guidance and manufacturer's declaration	64
[11.] accessories and options	66
[11thA] Vests for subjects for h/p/cosmos airwalk and h/p/cosmos airwalk se	66
[11.B] Wheelchair ramp	66
[11.C] h/p/cosmos running machines	66
[11.D] Safety multiple socket	66
[12.] Certificates	67
[13.] Disposal	67
[13.A] Disassemble and cut up	67
[13thB] h/p/cosmos unweighting devices	67
[14.] Contact.....	68
[14.A] Service department	68
[14.B] Sales department	68
[14.C] Headquarters	68
[15.] Appendix – Instruction and commissioning	69
[15.A] Calibrating the unweighting indications	69
[15.B] Instruction into general operation	69
[15.C] Instruction into functions	70
[15thD] Instruction in maintenance works and safety checks	72
[15.E] Advice and support documents	73
[15thF] Confirmation of commissioning and introduction / instruction	74
[15.G] Service report – cos15531	75
[15.H] Control protocol – cos15533	76
[15.I] Patients record	77

[15thJ] Protocol for maintenance on h/p/cosmos airwalk se..... 78

[3rd] Device description

Construction

The h/p/cosmos airwalk devices are made of powder coated steel. Two vertical sidebars carry the traverse with the deviating pulley. Like a crane the device is overlapping the treadmill in a height of approximately 2,5m.

The unweighting unit is integrated into the system. It may either be a mechanic (winch / spring), pneumatic (pneumatic cylinder) or electric (servo motor) unit performing the weight support. The unweighting force is transferred via a static rope to the unweighting bar. The patient is wearing a special vest, which is fixed to the unweighting bar.

In order to support the patient, the unweighting vest is similar to a climbing harness with a very high degree of comfort.

Application

In neurological rehabilitation it is important for the patient to start exercising as early as possible. Therefore an individual and optimal un-weighting system is crucial for the patient. The h/p/cosmos airwalk unweighting system supports a natural gait pattern. The single-point suspension allows dynamic vertical movement when walking and at the same time allows freedom in movement and body rotations where wanted. Additional fixation straps for further stabilization may be utilized if desired and if recommended for the patient. The unweighting, depending on the progress of therapy, can be adjusted via a hand unit.

[4.] **Safety, warnings, precautions, prohibitions**



[3.A] **Safety precautions, safety regulations, prohibitions and warnings**

- Anyone not involved in the training session must keep at least a safety distance of 2 meters away from any part of the system.
- The Levi unweighting unit must only be combined with the h/p/cosmos airwalk se. Don't use any unweighting device which has not been installed or checked by h/p/cosmos.
- Training must only be carried out with closed front door of the Levi unweighting unit housing.
- Check that the deflection pulleys are turning freely every time before using the equipment. Do not carry out training if the deflection pulleys are jammed. If a jam occurs during a training session, stop the session immediately.
- Training with the h/p/cosmos airwalk se may only be carried out by a trained user who has read the user manual. This person must be present at all times to supervise the training session when the system is being used for training.
- The unweighting device is approved for patients with a maximum body weight of kilograms. Heavier patients must not be trained with the h/p/cosmos airwalk se.
- Make sure, before opening the unweighting vest, that the patient is no longer being supported.
- Check the support bar, the unweighting vest and the rope for any damage before every training session. If one of these parts is damaged, treatment of the patient must not be started or continued.
- If there are faults or anomalies or if there is evidence of damage (e.g. unusual sounds), the training session must be stopped immediately and a h/p/cosmos authorized technician must be informed without delay.
- In the event of a power failure the patient must be lowered using the emergency release. Do not continue training until proper operation is guaranteed.
- If the patient can not be lowered for technical reasons (winch and emergency release are no longer working), do not attempt to remove the patient from the device on your own. Ask for assistance from one or more people you have called upon, and they can lift the patient while you then undo the buckles on the unweighting vest.
- The unweighting device may only be used if maintenance has been carried out properly according to the instructions included in this user manual.
- If a wheelchair ramp in combination with a treadmill (option) is used, make sure that the ramp is clean and correctly placed on the treadmill.
- You will find further safety instructions in part [6.M] of this user manual.

[3.B] **Danger zones**

No.	Illustration	Description
		<p>Illustration shows possible danger zones: there may be a risk of being jammed by incorrect and/or inattentive use.</p> <p>Please be careful and follow all safety instructions.</p>

[3.C] **Emergency release**

- To activate the emergency release, you must pull the knob out slowly.
- It is important to observe the patient all the time when doing this to ensure that he is lowered in a controlled manner. If possible, a second user should support the patient, while the other user is activating the release.
- Release the knob when the patient has been lowered to the desired level.






Diagram showing how to use emergency release


 **You must pull out the knob very carefully to enable you to control the speed at which the patient is being lowered (the speed when lowering heavy patients tends to be greater than when lowering light patients).**

 **A regular check (every three months) must be carried out to see that the emergency release is in good working order**

[3.D] **General Instructions / Safety**

This instruction manual is part of the device and must always be accessible for every user.

Exact observance of the instruction manual is a requirement for the appropriate operation of an h/p/cosmos device.

 **Important notes, warnings and safety precautions are marked with the opposite sign. It also reminds you of concerns which have to be considered for measurements and connection with other devices..**

The safety of your subjects and the compliance to the mentioned accuracy of measurement can only be achieved, if the consumption goods, described in the instruction manual, are being used and the recommended preventive maintenance and safety checks are performed by personnel authorised by h/p/cosmos.

Absorption-, cleansing- and disinfecting chemicals are only to be kept in the appropriate reservoirs. At the connection of gas pressure container it has to be examined about its MOT-admissibility (TÜV) first. Pull the plug before cleansing or disinfecting electric devices.

Safety, reliability and function can only be achieved if ...

- installation, extension, alteration and repair is performed by authorized people.
- the room for installation corresponds to the DIN and VDE installation directives.
- the device is plugged in at a socket with a protection system.
- the room for installation corresponds to the requirements of the surroundings for the device.
- the device after the instruction is being used in accordance with the instruction manual.

The entry of fluid into the device has to be removed immediately by the authorized customer service and a safety examination has to be performed. Damaged socket connections, wires and pressure control switches have to be replaced instantly by trained personnel or by authorised persons.



In case of any detected and/or assumed malfunctions and/or defects or unreadable safety warning labels the device has to be taken out of operation, the device has to be marked and secured against operation (for example pull power plug and affix warning/defect label on the power plug) and the supplier and authorised service personnel has to be informed in writing.

[3.E] Instructions for safety and operation

The regulations mentioned in this chapter refer to the Federal Republic of Germany. For all other states take into account the national legislation. If the running machine is a sports device (only **CE** mark without 4-digit number), it must not be connected with a medical product and must not be used for medical purposes. Due to higher safety requirements medical instruments and devices are not allowed to be connected via extension cable and/or multiple plug socket with the power supply, and a medical product must not be connected with a sports device either.



Disregard of warnings, disregard of intended and forbidden use, safety precautions and also unauthorized or lack of maintenance and/or regular safety checks may lead to injuries or even death and/or can damage the device and will result in loss of any liability and warranty

[3.E1] Electric safety / Safety classes

In order to protect the subject (athlete and/or patient) and the personnel, the association of German electro-technicians Inc, (VDE) has published special directives for medical used rooms and electro-medical devices. Devices with a power supply therefore have to, in order to prevent the passing on of the mains voltage over to touchable metal pieces, be equipped with not only a reliable isolation of the parts being under voltage but also additional safety precautions. The VDE-association divides it into so called safety classes.

Of the, for electro-medical devices, licensed safety classes are mainly used the safety class I (i.e. safety precautions with protective wiring), and the safety class II, (i.e. safety measures without protective wiring but with double isolation): Devices of the safety class I are devices, where the metallic casing-parts are connected with the protective wire of the line net via the safety contact. In the case of an isolation error the inserted fuse element switches off.

[3.E2] Patient security according to VDE 0750 (IEC 601-1) type B

If your h/p/cosmos device has a **CE** mark with a 4-digit number (e.g. **CE 0123**) at the identification plate, it is a medical product. At your h/p/cosmos medical product the patient is protected/isolated from the mains in accordance with the safety standard referring to the requirements for leakage current values according to IEC 601-1 type B. The connection of further mains-operated units to your h/p/cosmos-unit may cause that all the leakage currents add and the safety of the patient is reduced. Due to this the connection of further units may only be carried out on consultation with h/p/cosmos.

[3.E3] Potential equalization

The connection of e.g. EDP units to the V24 / RS 232 / USB / LAN / FireWire interface of the computer or of other independent measuring stations can only be carried out in medically utilized locations according to VDE 0107 / EN ISO 60601-1-1. Parallel connection of the potential equalization ensures, that in case of insulation failure none of the connected units exceeds a voltage of 10mV and thus endangers the patient.

For this purpose all of the h/p/cosmos units are equipped with a connection fixture with which the equipotent conductor, being part of the standard equipment, is connected to the equipotent bus bar of the room.

Standard potential equalization cable (length 5m) with 4mm² cross-section [cos10223]

[3.E4] **Connection of units installed outside medically utilized locations**

If units installed outside medically utilized locations (e. g. external printer, host computer etc.) are connected to an h/p/cosmos unit installed in a medically utilized location the VDE 0750 regulations must be observed. Connection only via ...

1. Optical fibre or optocoupler (4 kV checked and approved) or ...
2. Protective isolation via insulation transformer according to IEC 601-1, appendix k. These units are to be connected to the equipotent conductor.

Due to the increased safety requirements medical instruments may not be connected via extension leads or multi-way connectors.


[3.E5] **Medical electrical systems**

The medical electrical system (later in this text referred as "system") is a combination of several devices, of which at least one is a medical electrical equipment. A combination occurs when devices are connected via so called "functional connections" or via a mobile multiple socket. Not only electrical connections are defined as "functional connections" but also those for the transmission of signals, electrical energy and/or substances. Thus already the operation of two or more devices on the same multiple outlet is defined as a system. Even a mechanical connection or a wireless connection produces a medical electrical system.

Basically it is imperative that a system is as safe as a single medical electrical equipment. The requirements in the medical devices act and the applicable regulations (for example MPBetreibV) must be met. The system as a whole must ensure that within the patients environment or at contact with patients the same level of safety prevails, as determined by DIN EN 60601-1 (VDE 0750 Part 1).


[3.E6] **Subject's surroundings**

The surroundings of the subject must have a gap of 1.5 m as it has been proofed by experience. This has been laid down here as the surrounding.



[3.E7] **Protection against the danger of electrical shocks**

Casing and cover: Parts of non-medical electric devices within the surrounding of the subject, which, after having taken off the coverings etc. without the use of tools for a routine maintenance etc., are in touchable reach, have to operate with a voltage that does not exceed 25 volt alternate voltage and 60 volt direct voltage, which is produced by a separate source as described in IEC 601-1. According to this example the leakage current would flow from the electric device to the earthed subject via the trainer/doctor.



Electrical equipment [3] and a subject [1] must never be touched by the doctor or trainer [2] at the same time.

[3.E8] **Environmental requirements**

Unweighting devices are not to be used in medically utilized rooms with a danger of explosions or in easily inflammable atmospheres. The devices must not be installed near to e.g. an x-ray device, motors or transformer with high connection power, as the electric and magnetic interference can falsify measurements or even make them impossible. High voltage lines must be avoided nearby the device. h/p/cosmos electrical devices with mains connections must neither be used in wet and humid areas (e.g. swimming pool, sauna, etc.) nor in environmental chambers.

If not stated otherwise in the delivery information h/p/cosmos devices are designed for operation in normal climatic surroundings (DIN IEC 601-1):

- Temperature: + 10° ... + 40° C
- Relative humidity: 30 ... 70 % (non condensing!)
- Air pressure: < 1.5 bar (approx. 10,000 feet (3000m), without pressurization)

The system has to be protected from high humidity. Venting slots are not to be covered; otherwise it would hinder the air circulation. Store the devices at a temperature of – 20° ... + 50° C.












All devices can be stored without power-connection and without operation for a period of 6 – 9 months. After this period batteries in the device may be discharged and maintenance and/or a new setup/programming of the device could required by authorized h/p/cosmos technician.

[3.F] **Economic life-time**

The economic life time of the product is determined at common usage and application to 20 years, provided that after 10 years all electrical parts and component, as well as mechanical parts like for example deflection pulleys, are renewed (at damage or unusual heavy duty possibly earlier) and the recommended maintenance intervals are kept. Every maintenance and repair work needs to be carried out by authorized h/p/cosmos technicians. For expendable parts a shorter life-time is due.

[3.G] **Symbols used**

All symbols used comply with the according norms IEC417, IEC878, EN957-1:2005 and Council Directive 2002/96/EC.

	Read manufacturer's guide, advises, instructions and manual		Protection ground
	Warning / safety precautions Pay attention to accompanying instructions		Earth
	Devices of the type B		Potential compensation
	Dangerous electric voltage		Symbol for collection, treatment, recycling and disposal of waste electrical and electronic equipment (WEEE) as set out in Directive 2002/96/EC of 27 January 2003 of the European Parliament and of the Council on waste electrical and electronic equipment are necessary to reduce the waste management problems linked to the heavy metals concerned and the flame retardants concerned.
	Alternating current (AC)		
	Bullet point		
	Symbol based on ISO7010:2003-M002. Follow manufacturer's guide, advice, instructions in the manual. Manual contains relevant safety information.		

[3.H] **Safety standards, norms**

All h/p/cosmos unweighting devices are produced according to strict safety- and quality controls. The product must not be used, if the certificates/permissions and the (according to the list of technical specifications) defined safety standards do not correspond with the local and country-specific requirements. The country-specific requirements must be compared before initial operation. Operate unweighting device only, if the requirements correspond. Full list of applicable norms see also in the CE-certificate of compliance.

The standards with date in this manual with reference to design and construction characteristics, refer to the up to date valid versions at the time of manufacture.

The standards with date in this manual applying to routine measurements and continual safety relevant measures (e.g. electrical measurements to determine leakage current, insulation resistances etc.) can and will change with time. Therefore the up to date valid versions of the respective standards at the time of measurement and their limit values, which can deviate then from the procedures and limit values in this manual, apply.

[3.H1] **The CE mark**

The CE mark on the nameplate of the sports and fitness running machines confirms the compliance to the council directive 89/336/EEC appendix I (EMC electromagnetic compatibility). The audit has been conducted after the criterion of jamming and interference immunity.

[3.H2] **EMC Electromagnetic Compatibility**

The unweighting devices for medical applications (with CE 0123 mark) have been built corresponding to the European Union Council Directive 93/42/EEC

Norms: EMC: EN60601-1-2; EN61000-3-2; EN61000-3-3, Safety: EN60601-1

Certification bodies: RWTÜV Systems GmbH, D-45032 Essen

Strong electromagnetic fields, transmitters and interferences, which are above the normal tolerance, can interfere with measuring functions and displays of the running machines and can lead to other malfunctions.

[3.H3] **General advice**

The device is assigned to overvoltage-category II and pollution level I.

[3.] **Fields of application / intended use**

[3.11] **Professional application in sports and fitness**

The unweighting device h/p/cosmos airwalk se is a medical device. The user is free to use the system for sports and fitness applications.

See professional application in the medical field

[3.12] **Professional application in the medical field**

The purpose of using the h/p/cosmos airwalk se is to support treadmill training being used to treat patients with walking disabilities caused by cerebral, neurogenic, spinal, muscular or bone-related disorders.

The h/p/cosmos airwalk se can be used with the following patient conditions in particular:

Intended use / indications

Intended use

- Body weight support of a subject (during treadmill therapy/training)
- Fall protection of a subject (during treadmill therapy/training)
- Emergency stop in case of falling during treadmill therapy/training
- Balance training under unweighted and/or secured conditions
- Functional movement and gait training under unweighted and/or secured conditions
- Overspeed/hyperspeed and excess frequency training in athletics

Indications

According to: Johanna Jasper-Seeländer “Laufbandtherapie in der motorischen Rehabilitation”

- Hemiparesis
- Paraplegia
- Multiple Sclerosis
- Craniocerebral injury
- Further not inherent neurologic gait dysfunctions
- Parkinson's disease
- Rheumatism
- Orthopedically / surgically caused indications
- Amputations
- Complete Paraplegia with chip implantation
- Inherent gait dysfunctions

Just as for any other therapy, the physician in charge is always responsible for indicating treadmill training. The under section [3.14] listed contraindications must be taken into account for training on the h/p/cosmos airwalk se.

[3.13] **h/p/cosmos devices - MDD medical device directive. Safety standards**

Over 20 years of experience in the production of treadmill-ergometers and complex systems with accessories and special options have made h/p/cosmos the absolute treadmill specialist. We consider ourselves as the “market precursor in technology and safety” amongst all treadmill manufacturers. Safety has a very high priority for us –

because we know from experience what kind of accidents can occur and which consequences might result from these accidents.

An unweighting device is a very effective instrument in rehabilitation, therapy, diagnostics and other applications. However, like in many areas of daily life (i.e. all kinds of traffic, may it be by car, train or bus or the handling of electrical machine tools etc.) the use of a device such as an unweighting device involves dangers which can be minimized by abiding strictly to safety procedures, but can never be fully excluded. As a manufacturer it is h/p/cosmos' duty, and as the distributor or end-user/operator of the unweighting device it is your duty, to minimize those dangers as good as possible.

In medical applications with the combination of unweighting devices and running machines many additional risks have to be considered. These risks consider on hand to the connection to the system via the unweighting vest (fresh scars from surgeries, material in bones etc.) and on the other hand the therapy on the running machine and caused exercise (patients with pacemakers, ECG lead and electrode connections to the skin, health and physical limitations, high workloads during stress tests, medication, etc.). Therefore medical treadmills with interface and interactive heart rate based load controls are classified as risk class IIb medical devices. They have specific additional design characteristics compared to fitness treadmills, and can only be sold and serviced by authorized and certified distributors. h/p/cosmos has the responsibility to ensure that the systems are correctly installed (with usage instruction) and maintained. The unweighting device is classified in risk class IIa, but is in case of the conjunction to a running machine regarded as a IIb product.

Therefore we have implemented a medical device observation and reporting system based on the statutory obligations contained in, Guideline 93/42 EEC, MDD (medical device directive), ISO9001 and EN13485 standard and our internal safety policy and risk management system. This is of utmost importance and all our distributors must sign a contract that they will comply with these important international rules for medical devices, particularly for class IIb devices and their respective statutory requirements.

h/p/cosmos fulfils its duty in different ways:

- All medical unweighting devices have to undergo an internal product-type examination as well as external tests performed by the respective examination institutes such as the TUV Product Service in Munich / Germany or the EMC laboratory. Our medical treadmill-ergometers meet the norms and guidelines stipulated for medical applications such as the, MDD, EU guideline 93/42 EEC for medical products, IEC EN 60601-1 (VDE 0750), VDE 0751, EN 61000-6-1 and EN61000-6-3 (EMC) and EN 10535 (patient lifters)
- All h/p/cosmos products must undergo a detailed test after production. All tests are recorded, checked by a second independent employee and filed. h/p/cosmos treadmills have to pass the checks 100% without exception.
- h/p/cosmos has a quality management and risk management system according to ISO 9001:2000 & ISO 13485:2003 as well as medical device observation and incident reporting system.
- In the course of maintenance work h/p/cosmos technicians or authorised service personnel check the electrical safety of the unweighting device together with all safety relevant details and measure the accuracy of speed and elevation of the treadmills at the customer's premises.
- h/p/cosmos trains all relevant service technicians and sales personnel as medical device consultants with respective certification.
- The user manuals of h/p/cosmos unweighting devices contain all important information and safety advice as well as maintenance instructions and user regulations. The unweighting devices have a label with the date of the next safety check /maintenance.

Therefore we take it as our responsibility to inform our partners and customers about the background and necessity only to use medical products in the medical application field. Statutory requirements only permit us to supply medically classified equipment for medical applications.

Medical Device Directive Article 1 Definitions, scope 2 a.)

“medical device” means any instrument, apparatus, appliance, material or other article, whether used alone or in combination, including the software necessary for its proper application intended by the manufacturer to be used for human beings for the purpose of: Diagnosis, prevention, monitoring, treatment or alleviation of disease, diagnosis, monitoring, treatment, alleviation of or compensation for any injury or handicap, investigation, replacement or modification of the anatomy or of a physiological process, control of conception, and which does not achieve its principal intended action in or on the human body by pharmacological, immunological or metabolic means, but which may be assisted in its function by such means;

The use of a product not released for medical use by the manufacturer in a medical application can result in increased risks for the patients and for the operator – especially in the case of an accident on the treadmill ergometer. In this case all parties involved (your insurance company, your professional association, possibly patients lawyers and courts) will question negligence and liability. Therefore the use of a sports device for medical purposes bears increased risks for the distributor, customer, operator, end user and last but not least the patient. Despite that some of the above mentioned laws and guidelines are only binding within the European Union, the patient’s needs, the risks and the required safety procedures are the same all over the world. We ask all parties involved to help in improving safety.

[3.14] **Forbidden use - contraindications**

Like for any other therapy, the physician in charge is always responsible for indicating treadmill training. The following contraindications must be taken into account for training on the h/p/cosmos airwalk se:

Absolute contraindications:

- Acute myocardial infarction (within 2 d)
- Instable angina pectoris
- Cardiac arrhythmia pathology and/or limited hemodynamics
- Symptomatic massive aortic stenosis
- Uncompensated / uncontrolled heart insufficiency
- Acute pulmonary embolism or pulmonary infarction
- Acute endocarditis, myocarditis, pericarditis
- Acute aortic dissection
- Acute coronary syndrome
- Acute phlebothrombosis of the lower extremities
- Febrile infections
- Pregnancy
- Acute thrombosis
- Fresh wounds e.g. after surgery
- Acute fracture
- Damaged disc or traumatic disease of the spine
- Epilepsy
- Inflammations
- Acute migraine

Relative contraindications:

- Left main coronary stenosis
- Main artery disease
- Cardiac valve disease of moderate severity
- Known electrolyte imbalance
- Arterial hypertonia (RR > 200 mm Hg syst. > 110 mm Hg diast.)
- Tachyarrhythmia or bradyarrhythmia
- Hypertrophic cardiomyopathy and other forms of outflow tract obstruction
- Higher degree atrioventricular AV-blocking
- Anemia
- Physical and/or mental disabilities leading to inability to exercise adequately

Also see following guidelines:

Deutsche Gesellschaft für Kardiologie – Herz- und Kreislaufforschung e.V. (<http://leitlinien.dgk.org>).

American College of Cardiology Foundation - www.acc.org

American Heart Association - www.americanheart.org

http://my.americanheart.org/idc/groups/ahaecc-internal/@wcm/@sop/documents/downloadable/ucm_423807.pdf

Further contraindications may occur. This has to be evaluated by the responsible medical doctor.

The above list does not claim to be exhaustive. The decision as to whether a patient is suitable for treatment always comes under the remit of the physician in charge, who has sole medical responsibility for the treatment. As part of this, he must evaluate in particular, in each individual case, possible risks and side-effects of the treatment against the benefit gained from it. In addition, the patient's individual situation plays just as important a role as the basic risk assessment for specific patient groups.

Being a scientific discipline, medicine is subject to constant change in response to new knowledge and progress. It is therefore the task of the physician in charge to continually keep his knowledge up to date by reading the latest scientific literature and to acquire new knowledge during the course of treatment.



At the suspicion of unauthorized access or other reasons to lock the unweighting device, the system has to be locked, signed as „blocked“ and secured against use (for example, unplug the power cord and to affix warning label).

[5.] Installation & commissioning

[5.A] Transport, unpacking and packaging

When receiving the machine in a crate or unpacked, make sure the machine, the accessories and/or the packaging is not damaged. If you discover any damage and/or missing parts make a note on the packing-list / delivery note of the carrier. Inform h/p/cosmos and your dealer immediately in writing about any damages and/or missing parts.



The manufacturer does not undertake any liability for any damage, complaints and missing parts which is not reported immediately by delivery on the packing-list / delivery note.

Before you unpack the machine and accessories read instructions on the crate. Make sure that the machine, power connection cable or any optional equipment will not be damaged during unpacking. Pay special attention on small parts, so that you do not dispose them or any instructions with the packaging. Within Germany mostly all devices are delivered and assembled by h/p/cosmos directly or by an authorised forwarder. If delivered by h/p/cosmos the packaging will be taken back and recycled.

If the unweighting device is being delivered by a carrier, you can recycle the packaging yourself or send it back to the manufacturer (transportation is to be paid by the customer). Often a recyclable transport tool or packaging or a transport-fuse (metal angle with screws) is included in the extent of delivery.

Ask your dealer and the carrier to take the packaging and the recyclable transport tool back to your dealer or to h/p/cosmos at your own costs. In some cases a credit note can be granted. Special packaging and/or carrier constructions must not be disposed unauthorized.

[5thB] Transport to upper floors and through narrow doors

The unweighting device is usually delivered on a h/p/cosmos standard pallet (L: 275cm x W: 1,22cm x H: 170cm) with a crate. The pallet can either be brought in fully loaded through doors with a minimum width of 135cm with the help of a forklift truck or a handlift truck, or the single parts can be carried in by hand to the future location.



**Some of the components weigh between 45 and 80 kg.
Please notice, that the installation must be done by h/p/cosmos authorized technicians.**

Possible way of delivery on pallet in a crate.





Transports of heavy devices must only be carried out by authorised staff in compliance with the safety standards. Otherwise there is serious danger for people and devices.

[5.C] Mechanical installation

- In order to ensure proper installation and safety, either the manufacturer, an authorised service crew or an authorised dealer must always perform transport and installation of the devices
- For security reasons behind the unweighting device is a safety area of at least 2m x 1m length and width or a minimum of the width of the equipment must be free for space to move. Should a treadmill (option) be used in conjunction with a wheelchair ramp, in this case the safety zone is to measure from the end of the wheelchair ramp. In case of the option of reverse belt rotation, you must note that also the area in front of the treadmill is a safety area. So even in front of the treadmill, a safety area of at least 2 m x 1 m as a free fall space is required.
- The provided space for the unweighting device must be even and horizontal.
- Models with leveling sockets (adjustable "feet") at the rear of the unweighting device have to be adjusted so that they have a firm stand, otherwise it could lead to noises like knocking or rattling during the training. Check the level arrangement with a water-level at all axes on the frame.
- The bearing capacity of floor and ceiling in the building must be higher than the weight of the machine. It has to be approved for the h/p/cosmos device by an authorized body of the operator. Example h/p/cosmos airwalk see in standard configuration vor running machine h/p/cosmos mercury: Dimensions ground frame unweighting device: L: 2230 mm x B: 2070 mm = 4.62 m² platform. Net weight of unweighting device: 390 kg, additionally running machine: app. 200 kg; static bodyweight of subject: 200 kg, dynamical weight of subject: 1200 kg (up to 6 times the bodyweight), over all weight load of system at platform: 1790 kg = bearing capacity of platform: 387.45 kg (~ 400 kg) / m².

[5.D] Electrical installation



An overload or an voltage drop (even short-time) of more than 20 % of the mains voltage might cause malfunctions and/or defects and might totally switch-off the unweighting device.

- Installation of any electrical device and h/p/cosmos devices and running-machines must only be carried out at a voltage power connection with ground wire (earthing) including leakage-current protection-switch (interrupter) and according to VDE 0100 or/and the currently valid regulations and directives. Requirements for special locations, areas and establishments (e.g. medically used areas) must be strictly adhered. PE-connection (protection-earth-contact) is stipulated at all running-machines. The voltage drop between the beginning of the consumer's installation and the wall socket must not be more than 4% (DIN VDE 0100-520). It is the consignee's and user's personal responsibility to check the correct functions of the mains connection including the outlets and an authorized electrical engineer must check these points for perfect functions regularly (1 ... 4 years). Inspections of the electrical installations within the building are not incumbent on the supplier h/p/cosmos.
- As standard, a usual power supply of 230 Volt / AC, 50/60 Hz is sufficient for most of the running machines.
- Use a usual 16 Ampere expulsion fuse with C-tripping-circuit for your house distribution.
- For further questions please ask your electrical engineer or h/p/cosmos.
- Before installing the unweighting device please compare the specifications on the nameplate concerning the mains voltage and the mains frequency with your local characteristics. Connection only if identical.
- Check the main lead, the voltage power supply outlet and ground wire protection-contacts before plugging it in. Damaged leads and couplers and defective or dirty contacts have to be exchanged immediately. Rubber-leads can get porous and friable after some years.



- **Plug the device directly into the wall socket with checked ground wire. Each device has to be connected to a separate circuit. The socket has to be marked with name and serial number of the running machine. The use of an extension cables or a multiple plug sockets is not allowed.**
- **Electrical devices with mains connections must neither be used in wet and humid areas (e.g. swimming pool, sauna, etc.) nor in environmental chambers.**

[5.D1] Name plate



Example name plate h/p/cosmos airwalk se

For full details please refer to the chapter „technical data“ in this manual.

Devices produced until 01.01.2014 are labeled according with article number cos30017vaxx.

cos30017va01 was the airwalk se with treadmill support, cos30017va02 without treadmill support.

The new article number cos30017-01 is just the same as the former cos30017va02 (without treadmill support).

[5thD2] Electric safety measurements and „First Measured Values“

Immediately at first installation at the customer's site an electric safety test and measurement has to be performed for „Protective Earth Resistance“, „Isolation Resistance“ and „Leakage Current“ and the values have to be recorded on a special protocol [cos11690xx] and marked with „first measured values“. For details please refer to the chapter „maintenance / inspections“ in this manual. One copy of this protocol [cos11690xx] remains at the owner's manual and the original of the protocol with the „first measured values“ shall be sent to the manufacturer h/p/cosmos.

Picture: Example for electric safety testing device based on IEC 601-1



[5.D3] Potential equalization

The potential equalization cable [cos10223] must be connected with the connector on the device plug and with the potential-compensation-bearing within the medically used room. During installation, connecting or disconnecting the potential equalization, the device (unweighting device) must not be connected to the power supply. Protection against electric shock has to be provided in the finished product/device






- When used in the medical field, all devices within the system have to be connected via a potential-compensation-cable in „star-connection“ with the potential-compensation-bearing within the medically used room.
- First connect the potential compensation (potential equalization) with the corresponding plug pin (next to the main switch at the front) and then connect the mains plug.
- During electric safety measurements and test (leakage currents, etc.) the potential equalization cable has to be disconnected temporarily.


[5.E] Starting the unweighting device

[5.E1] Switching the device on

	If the unweighting device is being used with a running machine, each single switch on procedure has to be checked and coordinated to each other. If necessary priorities should be set.
---	---

No.	Illustration	Description
[01]		Plug in the power supply of the h/p/cosmos airwalk se unweighting device into the dedicated socket.

[5.E2] Switching the device off

	<ul style="list-style-type: none"> ■ The intervalls of switching Off- and On must not be shorter than 1-2 minutes. Otherwise it could lead to interference in the adjustment of the motor of the unweighting device or to the failure of the backup. ■ The system should be unplugged in case of longer breaks (for example over night) .
--	---

[5.F] Installation, Overview, Checklist

The authorized technician (distributor or service partner) has to check whether the installation of the device has been performed properly, before the customer is introduced into the operation of the h/p/cosmos device. For detailed information and instructions for installation and commissioning please refer to all respective individual chapters in this manual. Among others, following points save to be performed or controlled:















No.	Illustration	Discription
[01]		Compare the delivery note with the installed unweighting device. Are all parts like unweighting bar, rope, machine folder and accessories complete and without damage.







accomplished







Installation & commissioning






[02]		<p>After unpacking, lay components on a clean surface and screw in the levelling sockets according to the picture. screw in WITHOUT lock nut! DO NOT tighten the lock nut!</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
[03]		<p>After unpacking, lay components on a clean surface and screw in the levelling sockets according to the picture.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
[04]		<p>After unpacking, lay components on a clean surface and screw in the levelling sockets according to the picture.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
[05]		<p>After unpacking, lay components on a clean surface and screw in the levelling sockets according to the picture.</p> <p>Lay Components in right position for further assembling</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
[06]		<p>Connect parts with: 4 pcs. M8 x 30 8 pcs. NL8-Nord-Lock 4 pcs. nut M8</p> <div style="display: flex; justify-content: space-around;">   </div> <p style="text-align: right;">accomplished <input type="checkbox"/></p>




Installation & commissioning





[07]		<p>Turn around module for further assembling The module can now be placed on the future position of the system.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>						
[08]		<p>Lay Components in right position for further assembling</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>						
[09]		<p>Slide ground frame parts of the basic framework into each other and screw together.</p> <p>Screwing with 4 pcs. M8 x 16 panhead screws per side.</p> <p>If floor ist uneven, it can be compensated by the levelling sockets or the with the help of a screw driver (has to be removed afterwards).</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>						
[10]		<p>Slide frame parts for treadmill assembling into each other and adjust for the choosen h/p/cosmos running machine. see [11]</p> <p>Screw with 8 pcs. M8 x 16 panhead screws per side.</p> <p>If floor ist uneven, it can be compensated by the levelling sockets or the with the help of a screw driver (has to be removed afterwards).</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>						
[11]	<table border="0" style="width: 100%; text-align: center;"> <tr> <td style="width: 33%;">h/p/cosmos pulsar</td> <td style="width: 33%;">h/p/cosmos quasar</td> <td style="width: 33%;">h/p/cosmos mercury/locomotion</td> </tr> <tr> <td></td> <td></td> <td></td> </tr> </table>	h/p/cosmos pulsar	h/p/cosmos quasar	h/p/cosmos mercury/locomotion				
h/p/cosmos pulsar	h/p/cosmos quasar	h/p/cosmos mercury/locomotion						
								
[12]		<p>Screw in all levelling sockets as far as they contact the floor and the thread bolts don't endure the frames surface.</p> <p>Level the whole frame group with the levelling sockets:</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>						







<p>[13]</p>		<p>Straighten up left side arch, put it on the ground frame and screw it together with 8 pcs. M8 x 40 panhead screws.</p> <p>Until the arch is bolted up, an additional helper should secure it from falling</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[14]</p>		<p>Mount deflection pulley on right side arch: 4 pcs. M8 x 20 4 pcs. NL8-Nord-Lock</p>  <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[15]</p>		<p>Mount deflection pulley on traverse: 4 pcs. M8 x 20 4 pcs. NL8-Nord-Lock</p>  <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[16]</p>		<p>Straighten up right side arch, put it on the ground frame and screw it together with 8 pcs. M8 x 40 panhead screws.</p> <p>Until the arch is bolted up, an additional helper should secure it from falling</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>

<p>[17]</p>		<p>Mount traverse with :</p> <p>4 pcs. M16 x 35 4 pcs. NL16-Nord-Lock</p> <p>Tightening torque: 200Nm</p> <p>An additional worker ist needed.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[18]</p>	 	<p>Mount block for Levi-unweighting unit. 3 panhead screws M8 x 16</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[19]</p>		<p>Dismount the back board of Levi unweighting unit. Put it on block mounted under No. [18] and bolt it with: Ellen head screw 6 pcs. M8 x 16 included spring ring and lining disc.</p> <p>This operation must be done by two workers.</p> <p>Remount the back board after fixation of the Levi-Unit.</p> <p>Fix the potential equalisation cable on the back board!</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[20]</p>		<p>Connect the led outside potential isolation wire with the potential isolation bolt on the right side arch.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[21]</p>		<p>Thread the rope in pulleys and guidances of the unweighting device. This procedure ist identical with the initial installation and with later rope changings:</p> <ol style="list-style-type: none"> 1. Thread the rope through the top pulley: the rope must be inserted between the pulley and the marked bolt.

		<ol style="list-style-type: none"> Thread the rope through the middle pulley: the rope must be inserted between the pulley and the marked bolt. Feed the rope through the round opening in the Levi housing from the top. <p style="text-align: right;">accomplished <input type="checkbox"/></p>
[22]		<p>Wind the rope around the pulleys (① - ⑧) in the unweighting unit.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
[23]		<p>CAUTION: The rope must be inserted between the pulleys and the guiding bolts!</p> <div style="display: flex; justify-content: space-around;">   </div> <p style="text-align: right;">accomplished <input type="checkbox"/></p>

<p>[24]</p>		<p>Wind the rope around the rope drum (2 complete turns)</p> <p>Feed the rope through the rope fastening</p> <p>Fix the rope temporarily by slightly tightening both screws (you should still be able to move the rope)</p> <div style="border: 1px solid red; padding: 2px; display: inline-block; margin-bottom: 5px;">Rope drum</div> <div style="border: 1px solid red; padding: 2px; display: inline-block; margin-bottom: 5px;">Rope fastening with 2 screws</div> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[25]</p>		<p>Attach the support bar</p> <p>This step is not necessary within the first delivery, because it's premounted by the producer. For later services, rope changes and maintenance-/cleaning jobs release the rope from the support bar and feed it in again.</p> <div style="border: 1px solid red; padding: 2px; display: inline-block; margin-bottom: 5px;">Open zippers</div> <div style="border: 1px solid red; padding: 2px; display: inline-block; margin-bottom: 5px;">You can easily undo the rope thimble's locking device using your fingers and then remove the shaft.</div> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[26]</p>		<p>How to set correct height for the support bar:</p> <p>Push/pull the rope downwards through the rope fastening and measure the appropriate rope length as follows:</p> <p>Residual length with standard traverse: 121 cm</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>


<p>[27]</p>		<p>Fix and knot the rope:</p> <p>Fix the rope by tightening both screws firmly until the clamping plate makes full contact.</p> <p>Knot the rope as shown below and then tighten the knot.</p> <div style="border: 1px solid red; padding: 5px; display: inline-block; margin-top: 10px;">Rope fastening and clamping plate</div> <div style="border: 1px solid red; padding: 5px; display: inline-block; margin-top: 10px;">The knot must be as close as possible to the rope fastening.</div> <div style="text-align: right; margin-top: 10px;"> <p>Safety knot</p>  </div> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[28]</p>	 <div style="border: 1px solid red; padding: 10px; margin-top: 10px;"> <p>Raising a person with a body weight of approx. 80 kg: The support bar can be raised to within 2 cm of the pulley.</p> <p>CAUTION: Before releasing the support bar, lower the person down to an appropriate height.</p> </div>	<p>Check for correct rope guidance and carry out a functional test: Run through the Levi's entire range of motion with a load (person) attached using the remote control, monitoring the rope's movement, especially in the Levi housing.</p> <p>NOTE: Once the height of the support bar has been adjusted (task 10), a patient with a body weight of 80 kg can be raised until the support bar is approx. 2 cm below the pulley.</p> <p>Be aware that, if the patient lift is at its maximum position, the support bar can be pulled against the pulley with a maximum of 80 kg at lighter weights or if the dynamic range of motion is entered. This may damage the cushion of the support bar.</p> <p>For this reason, ensure that the support bar is not pulled against the pulley.</p> <p>In this case, correct the residual length as shown in [25] und [26]!</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[29]</p>		<p>Lock the housing of the Levi unit using the key.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>


<p>[30]</p>		<p>Screw locking nut (arrow) on thread of the levelling sockets with a socket wrench size 24.</p> <p>Put on the cap afterwards.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[31]</p>		<p>Check the wall socket: Optical check for damages, burned contacts, proper fixation and good ground contact (no paint, dirt, deformation)</p> <p>Check out the correct pin assignment of the wall socket with a socket-tester with direct indication.</p> <p>Check the electrical connection of the device: Direct connection to the wall socket with separate fuse, no danger of stumbling, no extension cable, no multiple plug socket, separate circuit for the treadmill. If a correct connection to the wall socket is not possible, note this on the delivery note, inform the customer and put the device out of operation (if necessary) and safe against start-up.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[32]</p>		<p>Immediately at first installation an electric safety test / measurement has to be performed for „Protective Earth Resistance“, „Isolation Resistance“ and „Leakage Current“ and the values have to be recorded on a special protocol [cos11690xx] and marked with „first measured values“. For details please refer to chapter „maintenance / inspections“. One copy of this protocol [cos11690xx] remains at the owner's manual and the original with the „first measured values“ shall be sent to h/p/cosmos. The potential equalization cable has to be connected AFTER the measurements.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[33]</p>		<p>In case of using the h/p/cosmos airwalk se with h/p/cosmos running machine the running machine has to be set on the ground frame of the unweighting device, the plates of the rear sockets have to be dismantled and set into the fixings of the ground frame, that fit to the treadmill model.</p> <div style="border: 1px solid black; padding: 5px; margin: 10px 0;">  <p style="margin: 0;">treadmills from other manufacturers: only if authorized by agency confirmation!</p> </div> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[34]</p>		<p>mounting the wheelchair ramp according to separate discription.</p> <p>For example h/p/cosmos locomotion med (option)</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>

[6.]Operation

[6thA] Controlements and displays: Manual control unit and housing unweighting unit.


[6thA1] Manual control unit of Levi unweighting unit

 Do not over stretch spiral cable. Tap buttons only slightly.
The keystroke is confirmed by a noticeable easing.



[6thA2] **Displays of the Levi unweighting unit and the connection for manual control**

The displays consist of 2 mechanical scales and a digital LCD display



[7.] Training

[7.A] General notes

[7.A1] Suitable clothing

Wearing appropriate clothing can reduce the patient's risk of skin irritation and lesions. We suggest patients wear tight-fitting pants made of soft cotton. Clothes made of thick, rough, fabrics with thick seams or very loose pants are unsuitable. Synthetic materials can cause mild burns on the skin due to friction.

Patients should wear flat, closed shoes on the treadmill. Where orthoses or, for instance, raised heels are required, suitable shoes must be worn to ensure a regular gait pattern.



**You can find further safety instructions in the front part of this manual.
Please always pay attention on these instructions.**



If users have long hair there is a danger of their hair getting caught in the treadmill during treatment. For this reason, we recommend that long hair is kept tied up at all times using a hairnet or other means.

[7.A2] Verwendete Laufbänder und Rollstuhlrampen:

Please always notice in case of using the h/p/cosmos airwalk se in conjunction with a treadmill the safety- and using instructions of the treadmill.

[7.A3] Before starting training

You must explain in detail to patients about what the Andago does before their first training session. Make patients aware of the benefits and risks of training.

[7.B] **h/p/cosmos airwalk Vests**

There will be supplied a unweighting vest with the h/p/cosmos airwalk se in standard configuration

[7.B1] **Putting on the Vest**

Note:

Before putting the h/p/cosmos airwalk vest on, make sure the patient is wearing appropriate clothing. Loose fitting and slick clothing will tend to allow vest to ride up.

When putting on the vest, remember that a tight fit provides greater support and feels more comfortable when unweighted. Avoid pressure under the arms and around the legs. The vest must support the person from the waist, not from the arms, ribs or hips. Always fasten the straps in the order given in these instructions. In the following steps the terms left and right refer to the left and right of the person wearing the vest. All straps should be fully extended and adjusted later to fit.



[7.B2] **Sizes**

Size	Waist Size	Vest Color	Order number
Small	56 to 81 cm (<22 to 32 in)	red	cos10095
Medium	81 to 112 cm (32 to 44 in)	blue	cos10096
Large	112 to 140 cm (44 to 55 in)	yellow	cos10097



[7.B3] Adjust Chest Straps

Put the left side of the cross over chest strap through the d-ring on the inside right of the belt, pull back and attach to the Velcro. Pull the right side of the cross over chest strap tightly over to the right, through the d-ring and press down securely on the Velcro. The shoulder straps will remain loose until connected to the unweighting bar.



The cross over chest strap is a positioning strap: it will fit differently on different body types. It should be tight enough to hold the chest in place while strapping the waist, but not tight enough to restrict breathing. You may need to make several adjustments to position the strap correctly for each person: with women, for example, the strap should be just below the breast because Unweighting will raise the vest up so that the four-way stretch fabric will protect and support the bust. Make sure that the back D-ring is at the center back in alignment with the spine.



[7.B4] Buckle Waist

Buckle up both buckles. Tighten up the top buckle by pulling first on one side of the adjusting strap and then the other, keeping the buckle centered. Tighten up the second buckle using the same method. You may have to reposition the cross over strap and the first buckle to get a snug fit. Do not tighten enough to restrict breathing. Make sure that the back D-ring is at the center back in alignment with the spine.



- **Note:**
- **The small (red) vest has only one waistband buckle, the medium (blue) and the large (yellow) have two.**



[7.B5] **Cuff Leg**

Leg straps should be all the way out. Each leg strap has a long front strap that attaches to the two adjustable buckles. On the vest fronts set the strap to the inside buckles for a small person and to the outside buckles for a larger person. The front strap should be centered from the waistband through the cuff, with the strap on the cuff in the center front of the thigh. The bottom of the leg cuff is set approximately 5 cm above the patella (knee). While holding the cuff in place (just above the knee), wrap the cuff and strap behind the leg then between the legs keeping the cuff taut. Pull on the strap to tighten comfortably. Press down on the Velcro to secure the strap. Repeat the procedure for the right leg.

[7.B6] **Adjust Leg Straps**

The back leg strap should be in a W position. Tighten the back leg straps by alternately pulling down on the adjustment straps attached to the side buckles. Do not tighten so much as to impede the normal leg action. The front leg straps are adjusted by pulling down on the buckle straps. It is not necessary to have the front straps snug. If the leg cuff rides up too much, reposition and retighten the cuff. Try to keep the front strap in the center of the thigh.



On the back of the vest are two elastic Velcro straps. These straps are used for gait control adjustments and do not need to be attached in order for the vest to be used for unweighting purposes.

[7thC] Attaching the Vest to the Unweighting Bar

Each shoulder strap at the vest has a D-ring that attaches with karabiners to the unweighting bar. Adjust the upright posture of the patient by moving the D-ring back or forth on the straps. With this individual adjustment the patient can be prevented of falling over to the front or the back. Open the Velcro strap underneath the D-ring and adjust it. by moving forth or back. Close the Velcro.



For more flexion, slide the D-ring backwards and fasten.

For more extension, slide the D-ring forward and fasten.

Before Unweighting anyone, be sure the waist belt is snug and that the safety snaps are closed completely.

The most important thing is, that the patient is feeling comfortable. Be sure that the vest is not wedged into the underarms and that it does not restrict breathing. Be sure that the leg straps are not tight enough to restrict circulation in the legs. Readjust the vest and single straps as necessary.

[7.C1] Correcting Vest Problems

If the vest rides up under the arms, reweight and pull the vest down. Have the person tighten the abdominal muscles, then retighten the waist straps and leg straps. If the chest strap is too tight during exercising, it has to be loosened. If the leg cuffs ride up, the cuff has to be fixed due to placing the front leg strap in the center of the thigh.

[7.D] Velcro-Gait Control Straps

The innovative gait control straps are found on the back of the h/p/cosmos airwalk vest. These elastic Velcro straps offer gait enhancement abilities never before obtained with traditional treatment. Helping to control internal/external foot rotation, stride length and hip rotation are just a few of the capabilities of the p/cosmos airwalk vest.

Some problems during the starting adjustments and some problems with the gait-belt are the internal or external foot rotation or an excessive stride length.

[7thD1] Internal Foot Rotation, Excessive Stride Length

The gait straps are attached to the plastic slides straight down and Velcroed to the back of the leg cuff. This helps to bring the foot into the proper position for correct stride.



[7.D2] External Foot rotation

The gait straps are Velcroed to the front of the vest across the buttock, onto the top of the opposite thigh and attached to the leg cuff.



This helps to correct the external rotation of the foot/hip. The gait control straps are also used in the rehabilitation of stroke and patients with Neuro muscular problems. They can be used to increase the stride length.

[7thE] Putting on the h/p/cosmos airwalk Vest for Wheelchair Patients

- Put the left side of the cross over chest strap through the d-ring on the inside right of the belt, pull back and attach to the Velcro. Pull the right side of the cross over chest strap tightly over to the right, through the d-ring and press down securely on the Velcro.
- While holding the cuff in place (just above the knee), wrap the cuff and strap behind the leg then between the legs keeping the cuff taut. Pull on the strap to tighten comfortably. Press down on the Velcro to secure the strap. Repeat the procedure for the right leg. The front strap on the cuff should be centered on the thigh.
- Buckle up both buckles. Tighten up the top buckle by pulling first on one side of the adjusting strap and then the other, keeping the buckle centered. Next tighten up the second buckle using the same method. You may have to reposition the cross over strap and the first buckle to get a snug tight fit. Do not tighten enough to restrict breathing. Make sure that the back D-ring is at the center back in alignment with the spine.

Reposition leg straps if necessary after waist has been strapped. The shoulder straps will remain loose until connected to the unweighting bar.

When strapping the waist area, fit the vest as tightly and as low as possible (have the person lift the shoulders if able.)



The patient is now ready to be Unweighted. Before Unweighting anyone, be sure the waist belt is snug and that the safety snaps are closed completely.

- Place the wheelchair under the h/p/cosmos airwalk. For treadmills, you may need to use a ramp to place the person under the system.
- Pull the bar down and connect the extension with the karabiners to the unweighting bar.
- Unweight the person from the chair and remove the chair. Adjust the extension as needed.


Most of the comfort problems encountered during unweighting are the results of improper adjustment of the vest.



[7thF] **Entering the unweighting device – supportign weight / standing up**



- Move the patient in his wheelchair onto the treadmill below the support bar.
- Check that all the buckles on the unweighting vest are fastened securely.
- Check the rope for possible damage and ensure that it is in a good condition.
- Hook the support bar onto the eyelets on the unweighting vest (see Fig. 9) and lift the patient out of the wheelchair into a standing position by keeping the relevant button on the Levi remote control unit pressed down (see p34)

	Before lifting up the patient the treadmill has to be in horizontal position!
	Make sure before lifting the patient that the support weight set on the Levi system is definitely less than the patient's weight!

[7thG] **Setting the support weight and range of motion**

	Make sure that the patient's wheelchair is located outside the 2-meter safety zone and that it cannot accidentally roll away.
---	--

[01]		<p>Set the support weight you want by pressing the relevant buttons on the Levi remote control unit.</p> <p>NOTE: If you do not know the patient's weight, first select a support weight that is definitely less than the patient's estimated weight. The patient can then be fully supported by pressing the relevant button on the Levi remote control unit (indicator in the  position). The weight can now be read from the patient scales and the support weight selected.</p>
------	--	---




[02]		<p>Now set the range of motion by pressing the relevant buttons on the Levi remote control unit. The patient is continually supported by the support weight previously selected only if the indicator lies in the gray area between 0 and //// during the entire walk movement siehe Abb [04].</p> <p>Now start the training session.</p>
[03]		<p>The support weight and range of motion can be adjusted at any time to suit the patient's requirements, without stopping the training session.</p> <p>You do this by pressing the relevant buttons on the Levi remote control unit (siehe Abb. [03] / [04])</p> <p>Note: If only small adjustments need to be made, just touch the buttons without keeping them pressed down.</p>

[7.H] General advice for training

The user should find these instructions helpful when providing treadmill therapy manually. Apart from drawing on the literature published on this subject, these instructions are also mainly based on the experiences of physiotherapists. This series of instructions is obviously not exhaustive.

When receiving manual treadmill therapy, patients should carry out as many activities on their own as is medically sensible, making the best possible use of their remaining functions. The following principle should be applied in order to ensure this happens: depending on the patient's training status, provide him with as little assistance and only support as much weight as is medically indicated. The purpose of the treatment is to achieve regular, symmetrical movement in both legs.

The following points in particular must be borne in mind during training:

	Ask patient how he is feeling in regular intervals Stop training session if the patient does not feel well.
	Have the straps on the unweighting vest become loose? Stop the training session and retighten any loose straps on the unweighting vest. Check to see whether the straps and buckles are twisted.
	Have the unweighting vest's buckles opened? Stop the training session and close the open buckles. Check that the buckles are still intact. If you do not detect any defects you can continue training. The training session must be ended if any buckles are defective.

[7.I] Unweighting / weight support

In order to achieve optimum activity from the patient's leg muscles during treadmill training, the level of body weight support provided should not exceed 30-40% of the patient's own body weight, provided this is not contraindicated by any other medical factors. In terms of the clinical factor you should use to determine the optimum body weight support, we recommend the point where the leg is sufficiently stabilized during the stance phase. You should bear in mind during the course of the training that the weight support must be reduced as quickly as possible according to the patient's training status.


[7.J] speed

The decision as to which speed should be used during the training session depends on the patient's abilities. In general, the increase in the speed of movements should be adapted individually to the progress the patient is making during treatment.

Selecting treadmill speeds that are deliberately slow may be beneficial to patients when they are learning or relearning to walk, especially at the start of treatment. In particular, being able to achieve a more effective and also more deliberate control of their movement can have a positive impact.

[7.K] **Elevation**

In order to provide sufficient variation as the treatment progresses, training can involve walking uphill (see picture).





	<p>Please notice:</p> <ul style="list-style-type: none"> ■ Start the training always in horizontal position of the treadmill ■ By changing the gradient angle you have to readjust the dynamic body weight support.
---	--

[7.L] **Trainingduration/ distance und training frequency**

Training duration and distance are determined by the patient's physical ability. We recommend short training sessions of 10-15 minutes to begin with, which can be adapted to the progress the patient is making during the course of treatment.

Current scientific knowledge about training stipulates as a requirement for maximum motor learning success that the organism has recovered as fully as possible (90-95% recovery) (Schnabel et al.(1997)). This is equivalent, for instance, to a rest period of 18 hours between two training sessions, which basically allows a training session every day using the h/p/cosmos airwalk se However, the decision about whether training can be carried out every day must always be made by the user, based on the requirements and abilities of the individual patient.

[7.M] **Training troubleshooting**

	<p>Is there a person or object inside the danger zone? Stop the training session immediately and make sure that there is no other person or object inside the danger zone before restarting.</p>
	<p>Do you notice the treadmill making jerky movements or unusual sounds? Stop the training session and carry out a functional check without a patient. Contact the treadmill producer.</p>
	<p>The patient can not be lifted out of his wheelchair? Check that the selected support weight is less than the patient's weight. If this is not the case, select a support weight that is definitely less than the patient's weight.</p>
	<p>The patient can not be lowered for technical reasons (winch and emergency release are no longer working)? Do not attempt to remove the patient from the device on your own. Ask for assistance from one or more people you have called upon, and they can lift the patient while you then undo the buckles on the unweighting vest.</p>

[7.N] After training

- Stop the treadmill.
- Push the wheelchair behind the patient.
- Lower the patient into the wheelchair by pressing the relevant button on the Levi remote control unit.
- Before undoing the unweighting vest, make sure that the indicator for the range of motion is in the position 0.
- Move the patient in the wheelchair off the treadmill via the ramp. It is essential that the user is standing behind the patient to slow down the wheelchair on the ramp.
- Now loosen the belt straps and remove the unweighting vest from the patient.
- Remember to enter the training data and update the patient record after every training session.

[7thO] At the end of a training day

- Switch of the treadmill according to the manufacturers orders.
- Pull the plug out of the socket, to separate the unweighting device h/p/cosmos airwalk se from the current.

[7.P] Patient documentation

You are recommended to create a patient record which you can use to note for every training session and patient the date, size of the unweighting vest, patient's weight, body weight support, the walking distance, speed, as well as any other comments and observations on each training session.

[7.Q] Patient becomes unconscious

A patient may have difficulties with circulation control when standing up from a seated position due to reduced muscle activity and paralysis of the vegetative nervous system. If the patient suddenly becomes unconscious the emergency physician must be informed and the patient laid down immediately or tilted backwards if sitting in a wheelchair. .


8. Maintenance and safety inspections

The h/p/cosmos authorised service engineers are happy to help you in the case of occurring problems.

A preventive maintenance can avoid problems in the future and is indispensable for the safety of such technical devices. Therefore ask for an annual preventive maintenance contract, which is highly recommended by the manufacturer, from our service department. Some basic regular maintenance and regular safety checks are obligatory!

A therapy system can only be safe if it is maintained on a regular basis and according to specified instructions

Maintenance must usually be carried out by h/p/cosmos or other by h/p/cosmos authorized persons every 12 months. Some preventive checks must also be carried out at shorter intervals.

	<ul style="list-style-type: none"> ■ You are not permitted to use the h/p/cosmos airwalk se if it is more than 12 months since the last maintenance inspection. ■ In case of any detected and/or assumed malfunctions and/or defects and/or unreadable safety warning labels the device has to be taken out of operation, the device has to be marked and secured against operation and the supplier and authorized service personnel has to be informed in writing. ■ Disregard of warnings, disregard of intended and forbidden use, safety precautions and also unauthorized or lack of maintenance and/or regular safety checks may lead to injuries or death and/or can damage the device and will result in loss of any liability and warranty. ■ Before intervening in the device for safety reasons switch the unweighting device off and pull the mains plug out. ■ During all maintenance works and safety tests make sure that no third parties are directly or indirectly in contact with the device under test and/or the technician performing the test. Keep a safety area of 2 m radius clear.
--	--

[8.A] **Preventiv maintenance**

Before switching the device on always check the circuit cable, plug, outlet socket and circuit entry of the device.

[8.B] **Immediate maintenance**

Immediate maintenance is necessary if:

- the device has been under high mechanical stress (push, power supply cable and/or interface cable defect through driving over it or pulling it)
- fluid has got into the device
- cable and/or connector plug have been damaged (therefore unplug cable from device and replace it)
- coverings and/or safety warnings have fallen off
- The rope or the deflection pulleys show sign of wear
- a defekt or malfunction of the device has been detected or is suspected

Only a properly and regularly serviced device is safe. The maintenance of the devices has to be performed by the service engineers authorized by h/p/cosmos, preferrably within the scope of a maintenance contract.

[8.C] **Regular inspections/examinations**

For use in sports and medical areas as well as for private, public and military applications refer to the date at the inspection sticker at your device.

To keep the condition of the device in due order, examinations have to be performed repeatedly. Follow manufacturers guidelines!

For the unweighting device h/p/cosmos airwalk se a maintenance interval respectively technical safety checks (STK) of one year (12 months) has been set.

These examinations are only to be performed by trained and authorized electric technicians.

The main inspection sticker on the device certifies also the inspection of the optional equipment and the accessories. However, inspection intervals for optional equipment and accessories (e.g. rope of h/p/cosmos airwalk unweighting device, chest belt for safety harness, compressors, etc.) can deviate significantly from inspection intervals of the main device. Read respective manuals for details.

For all required measurements and control steps consider also local requirements in your country and see detailed instruction and protocol form order no. [cos11690xx].

For maintaining the h/p/cosmos airwalk se there is a spezific detailed instruction and protocol form available.



[8thC1] Visual inspection dirt / damage – daily before training

- Carry out a general visual inspection to check for any damage to the treadmill and support system (particularly the rope). If any damage is suspected the h/p/cosmos airwalk se must not be used.
- Carry out a check for dirt, especially on the harnesses, leg straps, padding and Velcro fasteners and treadmill ramp (slipping hazard!). Dirty orthopedic components should be washed.
- Check the stitching and catches on the harnesses and leg straps. Defective orthopedic components must be replaced.

[8thC2] Emergency release unweighting device – every 3 months

- Functional test (enter in the “Maintenance activities report”).
- Lift up a weight or a person of a standard weight (e.g. 70 kg).
- Release the weight/person making several “interim” stops using the emergency release (see Part [3.C] of this manual). It must be possible to release the weight/person in a controlled manner. When the knob is released, the weight must be stopped again.

[8thC3] Running machine – follow manufacturers guidelines

- Maintenance according to the manufacturers guidelines.

[8thC4] **Changing the rope – min. every year**

The rope must be changed once a year (every 12 months). This is one of the most important safety measures of all. [spare part: cos100320]

Changing the rope is only allowed to h/p/cosmos authorized technicians!





Figure Levi unweighting unit

- [01] At the beginning of the rope – changing loosen the screws of the rope fastening and undo or cut off any knots at the end of the rope.
- [02] Pull rope out of the weight support pulleys.
- [03] Remove rope completely from the pulleys
- [04] Remove rope from support bar – see [4.E] – No. [24]
- [05] Afterwards thread the new rope into the unweighting unit according to [4.E] – startin at No. [20].

[8thC5] Spare parts


pos.	article description	units per machine	article code
1.	Levelling socket h/p/cosmos airwalk se	11	cos100313
2.	Hexagon nut M16 DIN936 flat nut	11	cos100616
3.	circlip NL8-Nord-Lock	4	cos100308
4.	Button head socket screw M8x16 A2	16	cos12634
5.	Button head socket screw, M8x40 V2A	16	cos11100
6.	pulley bottom side	1	cos100319
7.	pulley top side	1	cos100318
8.	mains cable 5m	1	cos100037va10
9.	Un-weighting rope - 7m - Ø8mm	1	cos100320
10.	snap hook with catch	2	cos11687-01
11.	discharge frame without hangers	1	cos100566
12.	Plastic cover plastic d=30mm, black	7	cos10535
13.	Cover cap plastic d=38mm black	12	cos10536
14.	cover cap plastik 60x60x3mm, black	8	cos100365
15.	Vest Medium f. h/p/cosmos airwalk (blue)	1	cos10096-01
16.	Vest Large f. h/p/cosmos airwalk (yellow) (optional)	1	cos10097-01
17.	Vest Small for h/p/cosmos airwalk (red) (optional)	1	cos10095-01
18.	Manual h/p/cosmos airwalk se EN	1	cos30017man-v1.3-en
19.	key for controll desk levi	1	cos100642

[8thC6] Protective Earth Resistance RPE measurement




- Resistance between housing and ProtectiveEarth-Connection.
- The low resistance pass is to be controlled according to VDE 0701/0702 (sports and fitness machines) or VDE 0751/IEC 601-1 (medical devices) by the protective-resistance measurement with a measuring device for the protective-earth-resistance measurement.
- The connecting lead is to be moved while measuring for at least 5 sec. If the resistance changes hereby it is highly probable that the cable or the connectors have been damaged. In this case the cable has to be replaced and the device has to be repaired immediately.

[8.C7] Isolation Resistance R_{iso} measurement



- Resistance between „bringed“ L+N and ProtectiveEarth-Connection.
- Make sure that all isolations, that are under stress of the mains voltage, are being seized. All switches and contactors should be connected.
- The measuring is to be performed with measuring devices for the isolation-resistance measurement according to VDE 0701/0702 (sports and fitness machines) or VDE 0751/IEC 601-1 (medical devices).

[8thC8] Equivalent (alternative) leakage current I_{EDL} measurement



- Impedance measurement, indicating the current in the Protective-Earth-Cable
- The measuring is to be performed by a measuring device for leakage current measuring according to VDE 0701/0702 (sports and fitness machines) or VDE 0751/IEC 601-1 (medical devices).
- The measurement is equivalent to the single fault earth leakage current based on IEC 601.

[8thC9] Electric Safety Tester an measurements



Picture: example for electric safety tester based on IEC601-1.

For all measurements read operation manuals of the measurement devices carefully and verify the values and intervals with the local guidelines and laws. An appropriate inspection record for all measurements and instructions are available at the manufacturers. See also detailed instruction and protocol form order no. [cos11690xx]. At first installation at the customer's site and before first commissioning of the unweighting device device, the „first measured values“ have to be determined and recorded on a special protocol [cos11690xx] and marked with „first measured values“.

During each preventive maintenance and/or regular maintenance and/or also after each repair work on the device (even if it was only mechanical repair work!) all above mentioned electric safety measurements and checks (see protocol [cos11690xx]) have to be performed again. All values have to be compared with the „first measured values“ and allowed range of tolerance specified according to VDE 0701/0702 (sports and fitness machines) or VDE 0751/IEC 601-1 (medical devices).

In case any or more of the measured values is/are not within the allowed range of tolerance specified according to VDE 0701/0702 (sports and fitness machines) or VDE 0751/IEC 601-1 (medical devices), the device has to be repaired until the values are within the specified allowed range of tolerance. In case the device can not be repaired accordingly and immediately, the device has to be marked and secured against operation and the supplier and authorized service personnel has to be informed in writing.

[8thC10] **Building installation: Electric Checks, Protective Earth Function, RCD, Leakage Current**

The customer/operator has the operator duty to check the building installation (electrical installations and stationary installations) in regular intervals according to regulation BGV A3 of the professional association respectively information of legal accident insurance every 4 years, on their own response and expenses, for functions and safety of the complete electrical building installation. At operating sites or rooms of special function - such as e.g. climate chamber, pressure chamber, installation with special endangering, etc. - the interval is prescribed one year or shorter if applicable.

In any case use ground wire plugs with tested ground wires only. Existent earth leakage current protection switches (so called RCD or "residual-current-operated protective device") at the building installation must be tested by the customer/operator monthly (button at RCD switch in fuse box of the building installation) for correct and proper function. As a rule the RCD switches at the fuse box are labelled with "press monthly".

Normally a button at the RCD switch has to be pressed for testing (not the "main handle"). With that an earth leakage current is simulated. The RCD switch must then switch-off the complete electricity supply. These tests must be carried out by the customer/operator at a time, when the normal operation is not disturbed, all electrical devices and computers are switched off and nobody is endangered.

Due to the test of the RCD switch it possibly comes, depending on the circuit of the building electricity, to a switch-off of the complete or a part of the power supply of the building.

In this connection it is important to stress, that these are not regulations, standards and safety checks specifically for h/p/cosmos unweighting devices, but current regulations for all electrical devices and building installations worldwide and applicable for all electrical devices with metal housing.

[9.] Trouble shooting



In case of any detected and/or assumed malfunctions and/or defects or unreadable safety warning labels the device has to be taken out of operation, the device has to be marked and secured against operation and the supplier and authorised service personnel has to be informed in writing.

[9.A] Mechanical / Noise problems

The h/p/cosmos airwalk se is an unweighting device with a spring mechanic. There may occur some noise during charging and discharging of the springs. If you perceive unusual noises:

- Check the mechanical device on overloading:
Have foreign objects (towels, cleaning cloths, etc.) come into the interior of the device?
- Check whether the noise comes from the treadmill.

If there is no improvement, please contact h/p/cosmos or your authorized h/p/cosmos dealer.

[9.B] Interference factor

[9.B1] Electrostatic Discharge

If the user moves around the devices they can be electrostatically charged with up to several thousand volts. If then the user touches a metal piece, keys or display, it can lead to an electrostatic discharge between the user and the device. Electrostatic discharges can in certain cases result in an interference of the device. Generally those electrostatic discharges are without harm for the user as well as for the device, but can be quite unpleasant. The main causes for electrostatic discharges are the choice of clothes, the sole of a shoe and the movement. Very dry air and many light fittings can also lead to the same results.

Solution: Try different clothes or shoes, humidify the air in the room, and switch part of the light fittings off. Please inform the manufacturer if you detect such interference.

[9.B2] Source of Interference

The devices should not be installed near to e.g. an x-ray device, motors or transformer with high connection power, as the electric and magnetic interference can falsify measurements. Very strong sources of interference (e.g. above the limit according to EMT) can influence the functions of the device.

High-tension power lines nearby and electrical devices without **CE**- sign and without a certificate of compliance for electro-magnetic-tolerance should be avoided as well.

Pay special attention to the EMC (electro magnetic compatibility) data and manufacturer's declaration in the chapter „technical data / EMC“ of this manual.

[9thC] Voltage at the device housing / electric shock

[9thC1] Isolation from mains supply

In order to isolate the device from the mains supply unplug the mains cable from the socket or from the device.

[9thC2] Open (interrupted) earth wire / ground – shown here on the esample of a running machine

In case the earth wire (ground) is not connected (for example due to defect in the building installation), on the metal frame of the device a voltage will be life via the „Y-noise-suppression-capacity“. This characteristic is very common with almost all electrical equipment with metal housing and EMC-fault clearance filters. The earth wire and/or contactor at the wall socket [1] in the building or at any other part of the AC-power line/feeder is open [2]. In this case a voltage of approximately 110 volt runs via the capacitor, of the EMC-fault clearance filter (to be found in the device) between housing and earth. Via the capacitors of the incorporated EMC-filters a voltage [5] of approx. 110 volts AC is live between frame [3] and grounding [4], e.g. via the floor or heating radiator system or other grounded parts or machines.



- If the device is flawless (at device category sports = without Potential-Isolation-Transformer) the current flow is significant noticeable in case of touching any bare metal parts (for example screws).
- If the device is flawless (at device category medical = with Potential-Isolation-Transformer) the current flow is almost not noticeable in case of touching any bare metal parts (for example screws).
- If the device is defective (at device category sports and also medical), e.g. because of a defective isolation inside the device, the current flow may be life-threatening high in case of touching any bare metal parts (for example screws).



Take the device out of operation immediately, unplug it and secure it against new start-up. Order an authorised electrician to repair the contactor circuit in the building or/and at the device.

Read and follow the directives and information concerning check of ground-wire terminals and functions in the chapters installation, maintenance and safety checks of this operation and service manual.


[9.D] **Electrical interference**


Interferences and error codes might be caused by problems with the power supply or because of not enough maintenance (e.g. pollution).

- Check voltage supply. Do not use extension cords or multi-way connectors. Connect the machine direct to the socket in the wall. Each machine should have an individual circuit.
- Check – if connected – the running machine
- Also contact problems at connectors (loose connections) caused by vibrations can result in malfunctions. So please check cables and connectors for loose connections.

[10.] **Technical specifications**

[10thA] **h/p/cosmos airwalk se**

	<p style="text-align: center;">h/p/cosmos airwalk se</p>
<p>Controls (Display & Keyboard)</p>	<p style="text-align: center;">1 LCD display, 2 mechanical indicators, 4 keys, mechanical emergency release</p>
<p>Power supply (standard, read nameplate)</p>	<p style="text-align: center;">230 V ~ f: 50/60 Hz In case of use in conjunction with a running machine, an additional socket with the requirements of the running machine is needed.</p>
<p>Power input</p>	<p style="text-align: center;">1840 VA</p>
<p>Overvoltage category</p>	<p style="text-align: center;">Category II: Life-Earth < 300 V transient over-voltage limit 2500 V</p>
<p>Fuse (standard, read nameplate)</p>	<p style="text-align: center;">10 A</p>
<p>Power consumption</p>	<p style="text-align: center;">max. 8 A</p>
<p>Elevation motor capacity</p>	<p style="text-align: center;">240 W</p>
<p>Safety Standards</p>	<p style="text-align: center;">CE 0123 IEC 60601-1, IEC 60601-1-1, IEC 60601-1-2, IEC 60601-1-6 DIN EN 12100-1, EN 957-1, DIN EN 62353, EN ISO 14971, ISO 9001, EN ISO 13485</p>
<p>Leakage current</p>	<p style="text-align: center;">0.2 mA</p>
<p>Isolation transformer</p>	<p style="text-align: center;">1 x 100 VA</p>
<p>Safety class / -categories</p>	<p style="text-align: center;">I / IP20 / B</p>
<p>Mode of operation acc. IEC 60601-1</p>	<p style="text-align: center;">Time of operation of 10% <i>2/18 = 2 minutes continuous use followed by 18 minutes not in use</i></p>
<p>Field of operation / Accuracy</p>	<p style="text-align: center;">sports S I and medicine, accuracy \pm 2 kg</p>
<p>Classification according to MDD</p>	<p style="text-align: center;">II a</p>

	Dimensions of the h/p/cosmos airwalk se in different combinations with running machines:					
	h/p/cosmos airwalk se without running machine	h/p/cosmos mercury med h/p/cosmos stratos med h/p/cosmos 150/50 locomotion DE med h/p/cosmos 150/50 locomotion force DE med	h/p/cosmos quasar med h/p/cosmos stellar med	h/p/cosmos pulsar	h/p/cosmos pulsar 3p	h/p/cosmos venus 200/75 (rs) if mounted in a pit
Required room height	2670 mm standard – traverse 2980 mm high traverse for patients up to 220 cm					
Dimensions of frame length x width x height	1960 x 2070 mm	2230 x 2070 mm	2430 x 2070 mm	2630 x 2070 mm	2630 x 2070 mm	2420 x 2070 mm plus control unit
Track access height from floor	...	240 mm	280 mm	280 mm	280 mm	...
Weight of device	275 kg	600 kg	780 kg	810 kg	850 kg	1075 kg
Weight load limits to floors (traffic load to DIN 1055 part 3)	App. values, depending on the equipment of the running machine					
	3.64 kN/m ²	3.90 kN/m ²	3.94 kN/m ²	3.69 kN/m ²	3.77 kN/m ²	4.58 kN/m ²
max. permissible user weight	In case of any further questions related to load limits onto floors, please ask our service department.					
max. user height	135 kg / 297 lbs					
Dynamic body weight support	200 cm at standard equipped system system height 267 cm (105.11") with an elevation of more than 10% use may be limited by the running-machine model, the size of the patient and the kind of (sportive) exercise. An optional higher traverse is available for patients up to 220 cm. system height with high traverse is app. 298 cm (117.32")					
Operation range	approx. 1...75 kg Due to manufacturing tolerances of the body weight support mechanism the maximum selectable support weight cannot exactly be named. However, the value of 75 kg can always be achieved..					
Support weight graduation	18 cm (vertical range of motion) at dynamic weight support.					
Unweighting indication (elektronik Display)	continuous					
Environmental conditions: Transport and storage	Resolution: 1/10 kg on patients scale/display of the actual weight support accuracy: app. ± 2 kg - Due to rope elasticity and friction of the rope pulleys (indirect measurement via pulleys) the weight measurement system has a relatively coarse hysteresis					
Environmental conditions: Operation	temperature: -20...+60 °C humidity: 20...90% - without condensation barometric pressure: 700...1060hPa					
combinable h/p/cosmos running machines	temperature: +10...+30 °C (deviation on request) humidity: 30...75% - without condensation barometric pressure: 700...1060hPa Maximum operating altitude: approx. 10,000 feet (3000m), without pressurization					
	h/p/cosmos mercury, h/p/cosmos locomotion, h/p/cosmos quasar, h/p/cosmos pulsar in their different version h/p/cosmos venus 200/75 in pit at ground level installation.					

Ask for further details and optional equipment or visit www.h-p-cosmos.com. Subject to technical alterations without prior notice. E & EO.

[10thB] **EMC electromagnetic compatibility. Guidance and manufacturer's declaration**

Table 201: Guidance and manufacturer's declaration - electromagnetic emissions (for all equipment and systems)		
The unweighting device is intended for use in the electromagnetic environment specified below. The customer or the user of the unweighting device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment — guidance
RF emissions CISPR11	Group 1 Class B	The unweighting device uses RF energy only for its internal function. Therefore, its RF emissions are low and are not very likely to cause any interference in nearby electronic equipment. The unweighting device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR11	Group 1 Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	


Table 202: Guidance and manufacturer's declaration - electromagnetic immunity (for all equipment and systems)			
The unweighting device is intended for use in the electromagnetic environment specified below. The customer or the user of the unweighting device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 0.5 cycle	<5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 0.5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the unweighting device requires continued Operation during power mains interruptions, it is recommended that the unweighting device be powered from an uninterruptible power supply or a battery. Caution! unweighting devices require high capacity UPS because of high capacity motor drive.
	40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles	40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles	
	70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles	70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles	
	<5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 5 sec	<5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 5 sec	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Not applicable	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE <i>UT</i> is the A.C. mains voltage prior to application of the test level.			

Please pay special attention also to the other EMC related chapters of this manual:

- a) Chapter: Safety, Installation, Operation, POLAR heart rate monitor: EMC
- b) Chapter: Error, Malfunction, possible interference

Table 204: Guidance and manufacturer's declaration -electromagnetic immunity (for all equipment and systems that are not life supporting)

The unweighting device is intended for use in the electromagnetic environment specified below. The customer or the user of the unweighting device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the unweighting device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	$d = 1.17 \sqrt{1/V * P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.17 \sqrt{m/V * P}$ for 80 MHz to 800 MHz $d = 2.33 \sqrt{m/V * P}$ for 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unweighting device is used exceeds the applicable RF compliance level above, the unweighting device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unweighting device.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 206: Recommended separation distances between portable and mobile RF communications equipment and the unweighting device (for all equipment and systems that are not life supporting)

The unweighting device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the unweighting device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unweighting device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (in Watt) W	Separation distance (in meters) according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.17 \sqrt{1/V * P}$	80 MHz to 800 MHz $d = 1.17 \sqrt{m/V * P}$	800 MHz to 2.5 GHz $d = 2.33 \sqrt{m/V * P}$
0.01 W	0.12 m	0.12 m	0.23 m
0.1 W	0.37 m	0.37 m	0.74 m
1 W	1.17 m	1.17 m	2.33 m
10 W	3.70 m	3.70 m	7.37 m
100 W	11.7 m	11.7 m	23.3 m

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

[11.] **accessories and options**

[11thA] **Vests for subjects for h/p/cosmos airwalk and h/p/cosmos airwalk se**

The vests offer high compatibility at best attainable freedom of body motion. They do not have undercutting straps during gait and have an optimum fitting due to the individually adjustable strap system.

4 different vests provide adaptation to all body sizes and forms. A special vest for children is also available. Additional straps for gait correction vest are included in the shipment. The unweighting is adjustable by the testing person itself.



h/p/cosmos airwalk vest red, waist 55-88 cm [cos10095]

h/p/cosmos airwalk vest blue, waist 81-112 cm [cos10096]

h/p/cosmos airwalk vest yellow, waist 112-145 cm [cos10097]

h/p/cosmos airwalk vest children [cos10112]

[11.B] **Wheelchair ramp**

When using the wheelchair access ramp pay particular attention to the gap between ramp and running belt that hands, hair, clothing or other items do not get caught. The ramp should not touch the running belt as this will lead to excessive wear. Ensure that the ramp is stable and cannot slip. Be aware of danger of slipping or tripping on ramp. Do not attempt to mount ramp from the side. Only original h/p/cosmos wheelchair access ramps should be used.



wheelchair ramp 150/50 [cos16186]

[11.C] **h/p/cosmos running machines**

See page 4 of this manual.

For further details contact your authorized h/p/cosmos dealer.

[11.D] **Safety multiple socket**

The stable safety multiple socket [cos100157] can be used for the safe connection of a medical h/p/cosmos treadmill with peripheral devices such as PCs, printer or monitors.

It meets all requirements for multiple sockets for creating connections of a medical electrical system within the patients environment.



safety multiple socket [cos100157]

As claimed in the norm EN60601-1 for medical electrical devices, the multiple socket is accessible only with the help of a tool. With this method an unintentionally made connection of further devices to the system is avoided, which lowers the risk of an exceedance of the maximum allowed patient-leakage-current of 0.5mA at the complete system.

[12.]Certificates

Certificates based on ISO 9001 and EN ISO 13485 and CE declaration of conformity see website:

<http://h-p-cosmos.com/en/safety/certificates.htm>

[13.]Disposal

By request and at the expense of the client h/p/cosmos might take over the disposal of old devices and devices not longer functioning. Please contact service@h-p-cosmos.com for a detailed offer. Note the information for possible disposal of the unweighting device parts or components through the client or a subcontractor.

The h/p/cosmos devices are marked with following sign / symbol on the name plate:

Symbol for collection, treatment, recycling and disposal of waste electrical and electronic equipment (WEEE) as set out in Directive 2002/96/EC of January 27, 2003 of the European Parliament and of the Council on waste electrical and electronic equipment are necessary to reduce the waste management problems linked to the heavy metals concerned and the flame retardants concerned.



h/p/cosmos EAR WEEE-Reg.-No. DE 42594388

[13.A] Disassemble and cut up

Use personal protective equipment, when cutting up material of any kind with the appropriate tools (eye-protection, dust-mask, etc.). Contact service@h-p-cosmos.com to receive the safety-data-sheet according to European Commission Directive 91/155/EEC for a material.

[13thB] h/p/cosmos unweighting devices

h/p/cosmos running unweighting devices consist among other things of powder coated and galvanized metals from different producers and qualities, stainless steel parts, aluminium parts, plastics, rubber, electronics with cables, boards and condensers as well as batteries. These materials can be recycled by handing them to the official municipal valuable substance collection or to authorized disposal partners of valuable substance disposal. Pay attention to the regulations of the disposal company.

[14.] **Contact**

For additional orders and technical enquiries please have the model type, the serial number and the date of installation of your unweighting device ready. In case of any further questions about delivery dates, service or maintenance, orders for consumables etc., please contact the corresponding phone, fax or email for qualified help.

[14.A] **Service department**

phone	+49 86 69 86 42 0
phone direct	+49 86 69 86 42 25
fax	+49 86 69 86 42 49
email	service@h-p-cosmos.com

[14.B] **Sales department**

phone	+49 86 69 86 42 0
fax	+49 86 69 86 42 49
email	sales@h-p-cosmos.com

[14.C] **Headquarters**

h/p/cosmos sports & medical gmbh
 Am Sportplatz 8
 DE 83365 Nussdorf-Traunstein
 Germany
 phone +49 18 05 16 76 67
 fax +49 18 05 16 76 69
 email@h-p-cosmos.com
 www.h-p-cosmos.com



Building 1 (top picture)
 h/p/cosmos development & production
 Am Sportplatz 8
 DE 83365 Nussdorf-Traunstein

Building 2 (picture below)
 h/p/cosmos sales & service
 Feldschneiderweg 5
 DE 83365 Nussdorf-Traunstein



[15.] **Appendix – Instruction and commissioning**


Once the installation of the unweighting device has been completed, the h/p/cosmos employee / h/p/cosmos partner starts with the instruction and commissioning of the device. It is important to include all people in the instruction and commissioning, who are going to work with the unweighting device. After the instruction is completed the instruction protocol is to be signed by the h/p/cosmos technician and all trained persons and sent back, together with the signed delivery note and the registration form, to h/p/cosmos.

h/p/cosmos zurück gesendet werden.




[15.A] **Calibrating the unweighting indications**


Maintenance on medical products and electrical devices may only be performed by h/p/cosmos authorised and trained personnel

We recommend a maintenance contract to secure the safe and proper use of your system.

<p>[01]</p>		<p>Calibrating the unweighting scales:</p> <p>A complete reset on „0“-kg is in cause of the construction of the device not possible. A rest weight of app. 0,3kg will always be displayed.</p> <p><i>Short-circuit terminal 10 (calibration) and terminal 15 (GND) with each other, using a piece of wire.</i></p> <p>→ The display will show the currently attached weight as 0 kg.</p> <p>Calibrating them with the support bar and unweighting vest attached: (note in the service list)</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
-------------	---	---





[15.B] **Instruction into general operation**

<p>I</p>		<p>Hand over directions for use.</p> <p>Inform about directions for use (always keep one within reach of all users).</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[01]</p>		<p> Refer to general warnings and safety precautions according to the directions for use. Place the print-out of the safety precautions (DIN A4 form in the delivery folder) close to the unweighting device.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>

[02]		<p>Especially point out the necessary safety-zone of min. L: 200 cm x W: 100 cm (or wider if running deck is wider) in case of fall.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
------	---	--




[15.C] Instruction into functions



[01]		<p>How to switch the unweighting device on and off.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
[02]		<p>How to set the support weight you want, by pressing the relevant buttons on the Levi remote control unit.</p> <p>NOTE: If you do not know the patient's weight, first select a support weight that is definitely less than the patient's estimated weight. The patient can then be fully supported by pressing the relevant button on the Levi remote control unit (indicator in the <u>////</u> position). The weight can now be read from the patient scales and the support weight selected.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>

<p>[03]</p>		<p>How to set the range of motion by pressing the relevant buttons on the Levi remote control unit. The patient is continually supported by the support weight previously selected only if the indicator lies in the gray area between 0 and  during the entire walk movement siehe Abb [04].</p> <p>Now start the training session.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[04]</p>		<p>Show how the support weight and range of motion can be adjusted at any time to suit the patient's requirements, without stopping the training session.</p> <p>You do this by pressing the relevant buttons on the Levi remote control unit (siehe Abb. [03] / [04])</p> <p>Note: If only small adjustments need to be made, just touch the buttons without keeping them pressed down.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[05]</p>		<p>Explain the Levi unweighting unit .</p> <p>Explain how to plug in the hand held remote control unit.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>


<p>[06]</p>		<p>Note the correct way to put on subject vest and connection to the support bar. See page 38 of this user-manuals.</p>  <p style="text-align: right;">accomplished <input type="checkbox"/></p>
-------------	---	--

[15thD] **Instruction in maintenance works and safety checks**

<p>[01]</p>		<p>Advise the customer to do regular inspection of the rope and pulleys. The manufacturer recommends maintenance contract by authorized h / p / cosmos service technician. Not or poor or from unauthorized personal done maintenance work and / or repair work or technical security controls mean danger for man and machine, can cause defects and lead to extinction and loss of any liability an warranty by the manufacturer.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[02]</p>		<p>Note to the Proper the adjustment of the levelling sockets. (in case of a later necessary adjustment)</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[03]</p>		<p>Note that the screw locking nut (arrow) has to be replaced afterwards on thread of the levelling sockets and tightened with a socket wrench size 24 again. Put on the cap afterwards again.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>

<p>[04]</p>		<p>Inform about the correct connection of the wheelchair ramp if used.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
<p>[05]</p>		<p>Inform about periodically recurring maintenance intervals and obligatory safety controls and the risks and dangers in case of disregarding these points.</p> <p>Offer and recommend maintenance contract through authorized and trained personnel.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>

[15.E] **Advice and support documents**

<p>[06]</p>		<p>Hand over and explain delivery papers (folder). Fill out registration form and send it to h/p/cosmos sales. Let the customer sign the delivery note and instruction protocol and send it back to h/p/cosmos sales immediately. Ask for brochure of customer and send it to the sales department at h/p/cosmos.</p> <p style="text-align: right;">accomplished <input type="checkbox"/></p>
-------------	--	--

[15thF] **Confirmation of commissioning and introduction / instruction**

By signing this protocol the authorised h/p/cosmos technician and the h/p/cosmos customer confirm the receipt and understanding of all warnings, safety precaution information, the performed instruction and commissioning according to form cos15228-01. The customer and user confirm the receipt of the listed devices including all accessories and options according to the h/p/cosmos delivery note. Disregard of warnings, disregard of intended and forbidden use, safety precautions and also unauthorized or lack of maintenance and/or regular safety checks may lead to injuries or even death and/or can damage the device and will result in loss of any liability and warranty. Please fill out the instruction protocol and send back to h/p/cosmos **via Fax to +49 / 86 69 / 86 42 49, via email to sales@h-p-cosmos.com or via post**

h/p/cosmos sports & medical gmbh Am Sportplatz 8 DE 83365 Nussdorf-Traunstein / Germany

h/p/cosmos sports & medical gmbh
Am Sportplatz 8
DE 83365 Nussdorf-Traunstein
Germany

customer's (end-user's) stamp / customer address:

h/p/cosmos device, model name	device serial number

instructor	name in clear block letters	h/p/cosmos dealer / technician	date and signature

instructed persons / customer / user / operator	name in clear block letters	position / function / department	date and user's signature

[15.6] Service report – cos15531



h/p/cosmos service report form no.:

<input type="checkbox"/> Repair	<input type="checkbox"/> Service	<input type="checkbox"/> Maintenance contract	<input type="checkbox"/> Warranty enquiry	<input type="checkbox"/> Good will bases
Company / organisation			Customer number	
Contact Person			Location, department, room #	
Street			Business hours	
ZIP code / city				
Phone / ext.			Device type	
Fax			Serial number (important!)	
email			Date of delivery	
www.			Distributor company name	

Detailed description of error and who detected the error. Without serial number of the device no further processing is possible.

Data of running machine or device, read out of MCU / User Terminal

OP 02: km	OP 03: h	OP 04: h	OP 05 (before update):
OP 34: Incr/10m	OP 35: km	OP 48: km/h	OP 99 (at MCU2/3) locked <input type="checkbox"/>

Following minimum actions have to be performed and checked. Confirm all performed actions with: X

Running machine cleansed inside / outside	Motor / fan grid cleansed	Light barrier cleaned
Driving belt tension / position adjusted	Light barrier fixation and adjustment	Speed + elevation sensor adjustment / fixation
Running belt oiled with ml	Running belt tension adjusted rotations	Running belt position adjusted: turns to
Elevation device tested	Elevation device greased (spindle + cogwheel)	All fixation washers of the roller axle checked
Function test mechanical and electrical	Function MCU / UserTerminal / Display	Test running (with person kg) min.
Pulse measurement <input type="checkbox"/> yes <input type="checkbox"/> no	Pulse control <input type="checkbox"/> yes <input type="checkbox"/> no	Function test pulse measurement
OP 47 interval display reset	Error memory frequency inverter <input type="checkbox"/> erased	Maximum speed: km/h

Notes:

Firmware update

The Firmware of the product updated to latest version _____ the latest version is already installed update rejected by customer

All spare parts needed: including order number as well as serial number for replaced and new parts

Amount	Spare part including date of production	Order number	Serial number replaced	Serial number new	Remark

Additional costs: none lubrication and cleaning material delivery costs spare parts fastening material stickers as following list:

cos10144 address sticker	cos10253-01 operation / safety instr.	cos10508 caution danger zone
cos10512 running belt adjustment	cos11668 safety warnings LT	cos11787 base sticker annual check
cos11880 unplug before opening	cos12976 warning triangle sticker	cos14543 top sticker year

Control check-up according to VDE 0701 / VDE 0751

The checks, necessary after every service (even mechanical service!) have been performed according to VDE 0701 (sports equipment), respectively VDE 0751 (medical device) after finishing the service with separate protocol form cos15533 & the following reported measurements:

Final report

The problem described above could be solved and the running machine / device can be used without limitations. Date for next check _____

The problem described above is not solved completely. Device can be used but following actions are required: _____

The device is defective and not safe! The device has been marked, put out of operation, secured against operation and must not be used at all!


Start working: o'clock	Working time: hours	<input type="checkbox"/> Distance flat rate for km
Driving distance (one way): km	Driving time (one way): hours	<input type="checkbox"/> Total distance <input type="checkbox"/> Distance proportional

Place, date	Stamp / Signature customer, name in capital letters	Signature technician, name in capital letters

Original to h/p/cosmos factory in Germany Pink print copy for technician Blue print copy for customer
 h/p/cosmos sports & medical gmbh Am Sportplatz 8 DE 83365 Nussdorf-Trautsain phone +49 / 86 69 / 86 42 0 fax +49 / 86 69 / 86 42 49 service@h-p-cosmos.com www.h-p-cosmos.com

[cos15531] service report form revision Feb. 09. 2006

[15.H] Control protocol – cos15533



h/p/cosmos control protocol for electric medical and sports devices

Customer- / device data

Control form for service report no.: _____
 Please fill in the service report - no if available:

Customer, ZIP-code, city _____

Device type _____ Serial no. _____ Construction year _____

General information

Manufacturer h/p/cosmos sports & medical gmbh – certified according to ISO 9001 and EN 46001, classification of devices DIN VDE I (one)

Manufacturer _____, classification of devices _____

Measuring instrument Gossen Metrawatt Secutest 0751 / 601 P HJS Eurotest 0701 – 0702 – 0751 other: _____

Measured main voltage _____ V Fuse breaker in the house installation: _____ A / Type: _____

Technician _____ Signature _____ Date _____

Notes: _____

Check

Separate device from the supply network. Separate connections to attachments (e.g. via RS232 interface). Remove data lines, earthing (potential equalisation). After measurements have been performed all connections have to be reconnected!

Check visually for damage: including all accessories. Device: engine compartment, mains connection lead with relief of pull-strain, ground wire connection and ground wire assemblies.

Check visually for damage: complete mechanical and wearing parts. Drive belt, belt spanner, running belt, elevation element with fixation screws, frame weld seams, fixation of screws and nuts, maintenance checked according to maintenance list.

Clean device and engine compartment. Remove stains and dust from cooling openings of the fan casing and the cooling ribs of the drive motor and the lift motor, as well as the louver and the perforated metal plate cover of the frequency inverter.

Ensure that the safety instructions are complete and displayed. Check the danger stickers, protection covers of the running machine, the motor cover and plastic lid of the installation channel in the engine compartment for their existence and for damage, replace if necessary. Micro-fuses, where applicable, which are accessible from the outside must be checked for the correct value and the correct labelling.

Check and adjust if necessary the belt re-entry zones in the back and at reverse belt rotation in the front as well. Gap less than < 8mm. According norm draft 60601-2-xx © IEC:200X 62D/479/NWIP 2003-05-18 and according to EN 957-1 (1999-02); see "test finger"

Latest version of the user manual for use by customers on site available. User manual version is compatible to the installed firmware at the running machine.

Users, patients and other third parties must keep a safe distance during measurements and must not touch the device.

	<input type="checkbox"/>	<input type="checkbox"/>	
Check according to DIN VDE 0701, Edition 2000-09 All h/p/cosmos devices of the category sports – CE As well as all medical devices supplied before 14.06.1998 (MDD becomes effective). Check according to DIN VDE 0751 Edition 2001-10 All h/p/cosmos devices of the category medical – CE 0123	Limit values DIN VDE 0701	Limit values DIN VDE 0751	Result of measurement
Protective-resistance for devices with solid mains connection lead Check: device with mains connection lead	< 0.3 Ω	< 0.3 Ω	_____ Ω
Protective-resistance for devices with removable mains connection lead Check 1: measurement only at mains connection lead.	< 0.1 Ω	< 0.1 Ω	_____ Ω
Protective-resistance for devices with removable mains connection lead Check 2: measurement only at device.	< 0.2 Ω	< 0.2 Ω	_____ Ω
Isolation-resistance Measurement at U _{ISO} > 500 Volt	> 1.0 MΩ	> 2.0 MΩ	_____ MΩ
Isolation-resistance (medical devices supplied before 14.06.1998) Measurement at U _{ISO} > 500 Volt	> 2.0 MΩ	-----	_____ MΩ
Alternative device leakage current Test voltage 250 V	< 3.5 mA	< 1.0 mA	_____ mA
Earth- / device leakage current I_{SL} (As alternative to alternative device leakage current) Test voltage 250 V	< 3.5 mA	< 0.5 mA	_____ mA
Running machine / device is plugged into the wall socket directly. The use of extension cables or multiple plug sockets is not permitted. Running machine / device is installed according to safety regulations, a clear safety zone of min. L 2m x W 1m behind the device exists. Function check: accessories & device. Start device, elevate lift device, check all existing emergency-stops, test pull rope & safety-stop. Use also service report form cos15531. Check accessories according to separate maintenance list of options & accessories.			
Assessment of the check: checks and measurement results OK. Place test badge with date-code for next annual check.			
Assessment of the check: checks and measurement results not OK. Customer informed and device put out of operation and secured against use. If the device was put out of operation this has to be reported on the service report and must be signed by customer.			

[cos15533] control protocol revision Feb. 09, 2006

Original to h/p/cosmos factory in Germany Pink print copy for technician Blue print copy for customer
 h/p/cosmos sports & medical gmbh Am Sportplatz 8 DE 83365 Nussdorf-Trautstein phone +49 / 86 69 / 86 42 0 fax +49 / 86 69 / 86 42 49 service@h-p-cosmos.com www.h-p-cosmos.com

[15thJ] Protocol for maintenance on h/p/cosmos airwalk se

Maintenance activities report	page:
--------------------------------------	-------------

Device installed:	Date:	Company /Distributor):.....
Serial-No. h/p/cosmos airwalk se:	Name (Technician):	

Date	Maintenance activities (based on page 51 section [7.C1]visual inspection section [7.C2] emergency release	Defects	Action to be taken (e.g.: rectify defect / notify h/p/cosmos distributor / order replacement)	Stamp / Signature
		<input type="checkbox"/> no <input type="checkbox"/> yes		
		<input type="checkbox"/> no <input type="checkbox"/> yes		
		<input type="checkbox"/> no <input type="checkbox"/> yes		
		<input type="checkbox"/> no <input type="checkbox"/> yes		
		<input type="checkbox"/> no <input type="checkbox"/> yes		
		<input type="checkbox"/> no <input type="checkbox"/> yes		
		<input type="checkbox"/> no <input type="checkbox"/> yes		
		<input type="checkbox"/> no <input type="checkbox"/> yes		
		<input type="checkbox"/> no <input type="checkbox"/> yes		
		<input type="checkbox"/> no <input type="checkbox"/> yes		



[15thK] Connection Diagram

